

Genome editing made easy: Top 10 tips & tricks



Improve your genomic cleavage detection (GCD) assay

- 1 Optimize cell culture conditions**
 - Pay attention to passage number, seeding density, and media.
 - Find the best culture reagents for your cell line at [thermofisher.com/cellculture](https://www.thermofisher.com/cellculture).
- 2 Optimize transfection conditions**
 - High transfection efficiency and minimal cellular toxicity will maximize efficient genome modification and improve genomic cleavage detection analysis.
 - Find the best transfection reagent for your cell line and optimize transfection conditions with our [transfection reagent selection guide](#).
- 3 Optimize the amount of nuclease**
 - Find the best quantity of nuclease for your cell type and culture conditions for [TALs](#) and for [CRISPRs](#) in the corresponding product user manuals.
- 4 Include negative controls**
 - Plan to include a mock transfection with an irrelevant nuclease as a negative control to distinguish background from cleavage banding in the gel analysis portion of the [GeneArt™ genomic cleavage detection assay](#).

- 5 Prepare cell lysate**
 - Ensure that the cell lysate contains at least 50,000 and no more than 2×10^6 cells in a 50 μ L cell lysis reaction—lysate from less than 50,000 cells will be insufficient for PCR, and lysate from more than 2×10^6 cells will inhibit PCR.
- 6 Generate a strong PCR product**
 - Generate a reaction that produces a single PCR product that migrates on a gel as a crisp band with minimal background.
 - We recommend designing a PCR fragment approximately 500 bp in length to allow for analysis of cleaved products.
 - Include a PCR sample from the mock transfection with the irrelevant nuclease (see Tip 4).
- 7 Verify PCR product before proceeding**
 - A single band is crucial for obtaining accurate cleavage detection. If a single band of correct size is present with intensity similar to 50 ng of 400 bp band in the mass ladder, proceed to the denaturing and reannealing step.
- 8 Denature and reanneal PCR fragment**
 - It is also good practice to include a negative control sample that will not be treated with the GCD enzyme, to distinguish background from true cleavage banding.
- 9 Optimize enzyme digestion**
 - To prevent over-digestion of the sample, be sure not to exceed the incubation time with the GCD enzyme, and continue immediately to gel loading.
- 10 Optimize gel analysis**
 - Optimize the gel imaging conditions to achieve the best signal/background ratio.

Find out more at thermofisher.com/genedetect

ThermoFisher
SCIENTIFIC