

# In-depth Gene Expression Profiling of Seminomatous Testicular Germ Cell Tumors

Sunali Patel<sup>1</sup>, Shirley Chu<sup>1</sup>, Lovorka Degoricija<sup>1</sup>, Kathy Y. Lee<sup>1</sup>, Ad J.M. Gillis<sup>2</sup>, Martin Rijlaarsdam<sup>2</sup>, Lambert C.J. Dorssers<sup>2</sup>, Leendert H.J. Looijenga<sup>2</sup>  
<sup>1</sup>Thermo Fisher Scientific, 180 Oyster Point Blvd, South San Francisco, CA, 94080, USA and <sup>2</sup>Erasmus Medical Center—University Medical Center Rotterdam, P.O. Box 2040, 3000 CA Rotterdam, The Netherlands

## ABSTRACT

In depth Whole Transcriptome Analysis by RNA-Seq was used to compare the expression profile of testicular germ cell tumors of seminoma histology and normal testis (Poster #2, Abstract #3561). From this analysis we identified a large number of genes that showed differential expression with statistically high significance, ( $p < 0.01$ ,  $> 2X$  FC). The list included genes related to stem cell pluripotency (*NANOG*, *POU5F1*), proliferation (*KRAS*, *CCND2*), a number of non-coding RNAs (snoRNAs, miRNA precursors, *XIST*), and testicular cancer related genes (*LDHB*, *AKAP4*). From this list we identified a subset of genes for validation and further screening. We used the QuantStudio™ 12K Flex Real-Time PCR OpenArray® platform to quantitatively screen up to 112 coding and non-coding genes (including controls) using samples from different malignant testicular germ cell tumors (GCC): seminoma (SE) and non-seminoma (NS) and normal tissue (N).


We used a complete work flow solution from sample prep to next-generation sequencing (NGS) to qPCR to compare the expression profile of normal testis and seminoma type germ cell tumors.

## MATERIALS AND METHODS


### A. Profiling on NGS Ion Proton™ Platform



### B. Validation/Screening on OpenArray® Platform



**Figure 1. Complete Workflow.** The complete workflow from profiling with sequencing on the Ion Torrent™ Proton™ System (A) to screening and qPCR validation on the QuantStudio™ OpenArray® Platform (B).




**Figure 2. The QuantStudio™ 12K Flex Real-time PCR OpenArray® System.** Custom 112 format OpenArray® plates are used to screen 112 genes for each sample. 24 Samples (14 SE, 7 NS, and 3 N) were screened in triplicate on four custom OpenArray® plates. 200ng of total RNA (converted to cDNA) was used to allow ~1ng per through-hole qPCR reaction.




**Figure 3. Distribution of 112 genes spotted on QuantStudio™ 12K Flex Real-time PCR OpenArray® Panel.** Genes selected had consistent differential expression between normal and GCC samples, and are suggested to play a role in normal and aberrant gonadal development.


## RESULTS



**Figure 4. Cluster Analysis of OpenArray® Panel Screening Data.** Expression levels were normalized to USP32 (most stable in samples according to transcriptome sequencing and qPCR data), and hierarchically clustered using an unsupervised centroid algorithm. Green means  $\Delta\Delta Ct$  is greater than the median of that gene in all samples; red means  $\Delta\Delta Ct$  is lower than median; grey is undetectable. SE samples are separated from NS and N samples using these selected genes.



**Figure 5. Comparison of Results using Open Array and NGS.** Fold change comparison between normal (3 samples) and SE (3 samples) that were run on Open Array and NGS for 112 genes. qPCR expression levels were normalized with USP32. 94% of the assays showed good directional concordance. Blue – OpenArray® panel FC (fold change) for qPCR; Red, TopHat, FC for NGS.



**Figure 6. Comparison of FC Across Three Different Platforms.** A subset of 20 TaqMan® Gene Expression Assays comparison between 3 platforms: NGS (STAR and TopHat2/BowTie2 aligners), TaqMan® Assays on 384-well plate and OpenArray® panel. qPCR expression levels were normalized to USP32, HPRT, and RNPS1.

## CONCLUSIONS

1. We successfully used the Ion Proton™ System with Ambion® RNA library preparation kits for Whole Transcriptome sequencing on normal and seminoma cancer samples (Poster # 2, Abstract 3561) following by gene expression validation and screening on QuantStudio™ OpenArray® platform (Figure 1 and 2).
2. From NGS data we identified 112 genes of interest that had  $p < 0.05$  and  $> 2$  fold differential for screening on the QuantStudio™ 12K Flex Real-Time qPCR Systems OpenArray® Platform (Figure 3).
3. Data from 24 samples show that the non-seminoma, seminoma, and normal samples cluster separately and can be differentiated (Figure 4).
4. Fold change across three different platforms, NGS and OpenArray® and TaqMan® Assays on plates showed good directional concordance (Figure 5 and 6).

## ACKNOWLEDGEMENTS

We would like to thank Dr. Stephen Jackson for his great support with data analysis and generating the heat map.

## TRADEMARKS/LICENSING

The QuantStudio™ 12K Flex Real-Time PCR System is For Research Use Only. Not for use in diagnostic procedures. © 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. TaqMan is a registered trademark of Roche Molecular Systems, Inc., used under permission and license. Partek is a registered trademark of Partek, Inc.