

Discover the connection

ImmunoCAP™ Allergen Components testing

Whole Allergens and Allergen Components help you diagnose allergy, allowing you to prepare a more comprehensive management plan.

Optimize management to help:

- ASSESS risk for systemic allergic reactions
- EVALUATE potential reaction to baked egg or milk products
- ADDRESS parental anxiety

With **ImmunoCAP Allergen Component test** results, you have more of the information necessary for proper diagnosis, allowing you to evaluate your patient's potential risk of systemic reaction, and develop a more comprehensive management plan.



PEANUT



Peanut Allergen Component testing can help determine which proteins your patient is sensitized to.

A specific IgE blood test that detects sensitization to the whole peanut is the first step in discovering the likelihood of a systemic reaction and the necessary precautions that may be prescribed.



77.6%

The percentage of patients sensitized to peanut who may not be at risk for a systemic reaction.¹

Characteristics of individual proteins

Peanut f 13	<ul style="list-style-type: none"> • High levels of peanut IgE can predict the likelihood of peanut sensitivity, but may not be solely predictive of reactions or allergic response¹
CCD MUXF3	<ul style="list-style-type: none"> • LOWEST RISK of systemic reaction² • Highly cross-reactive with pollen, plant food and venoms²
Profilin Bet v2	<ul style="list-style-type: none"> • LOWER RISK of systemic reaction^{2,3} • Cross-reactive with pollens²
Ara h 8 f 352	<ul style="list-style-type: none"> • LOWER RISK of systemic reaction^{4,5} • Risk of mild, localized symptoms, such as itching/tingling of the lips, mouth, and oropharynx⁴ • Cross-reactive with pollens (e.g., birch)⁵
Ara h 9 f 427	<ul style="list-style-type: none"> • VARIABLE RISK of systemic reaction including anaphylaxis^{7,8} • Often accompanied by sensitization to other peanut proteins⁹ • Cross-reactive with fruits with pits (e.g., peaches)⁷
Ara h 1, 2, 3, 6 f 422, f 423, f 424, f 447	<ul style="list-style-type: none"> • HIGHER RISK of systemic reaction including anaphylaxis^{7,10} • Sensitization to Ara h 2 is nearly always associated with clinical peanut allergy¹¹

CCD, Profilin, Ara h 8 MUXF3, Bet v2, f 352	Ara h 9 f 427	Ara h 1, 2, 3, 6 f 422, f 423, f 424, f 447	Management Considerations ^{1,4,12-18}
+	-	-	Oral food challenge (OFC) with a specialist may be recommended. High likelihood that patient may pass OFC. If patient passes an OFC: <ul style="list-style-type: none"> • Foods prepared with or around peanuts may be consumed • Patient not restricted to peanut-free zones
+/-	+	-	<ul style="list-style-type: none"> • If there is no clinical history of symptoms, please see considerations above • If there is a clinical history of symptoms, please see considerations below
+/-	+/-	+	<ul style="list-style-type: none"> • Choose peanut-free zones for patient's safety • Consider prescribing epinephrine auto-injector • Family, colleagues, and teachers should be made aware of allergy and have a plan

As in all diagnostic testing, any diagnosis or treatment plan must be made by the clinician based on test results, individual patient history, the clinician's knowledge of the patient, as well as their clinical judgment.

EGG



Egg Allergen Component testing can help determine which proteins your patient is sensitized to.

A specific IgE blood test that detects sensitization to egg white is the first step in discovering your patient's allergy. Egg Allergen Component tests can help you determine the likelihood of reaction to products baked with egg, such as muffins or cookies, as well as the likelihood of allergy persistence.



70%

The percentage of children with egg allergy who do not react to baked egg.¹⁹

Characteristics of individual proteins

Protein	Gal d 2 / f 232	Gal d 1 / f 233	Management Considerations ^{19,21}
Egg White f 1			<ul style="list-style-type: none"> Avoid uncooked eggs Likely to tolerate baked egg Baked egg oral food challenge with a specialist may be appropriate Consider repeating IgE component test biennially during childhood to determine potential tolerance May be transferred via breast milk, so mothers of infants with egg allergy should take caution when breast-feeding
Ovalbumin	+	-	
Ovomucoid	+/-	+	

*In clinical studies, extensively baked muffin and waffle were heated to the point of protein denaturation.

Protein	Gal d 2 / f 232	Gal d 1 / f 233	Management Considerations ^{19,21}
Ovalbumin	+	-	<ul style="list-style-type: none"> Avoid all forms of egg Consider repeating IgE component test biennially during childhood to determine potential tolerance Patients sensitized to ovalbumin with low levels of IgE to ovomucoid may react to egg that is not fully baked
Ovomucoid	+/-	+	

As in all diagnostic testing, any diagnosis or treatment plan must be made by the clinician based on test results, individual patient history, the clinician's knowledge of the patient, as well as their clinical judgment.

MILK



Milk Allergen Component testing can help determine which proteins your patient is sensitized to.

A specific IgE blood test that detects sensitization to cow's milk is the first step in discovering your patient's allergy. Milk Allergen Component tests can help you determine the likelihood of reaction to baked goods, such as cookies or cheese pizza, as well as the likelihood of allergy persistence.



75%

The percentage of children with milk allergy who do not react to baked milk.²⁸

Characteristics of individual proteins

Protein	Bos d 4 / f 76	Bos d 5 / f 77	Bos d 8 / f 78	Management Considerations ^{28,30-38}
Cow's milk f 2				<ul style="list-style-type: none"> Avoid fresh milk Likely to tolerate baked milk products Baked milk oral food challenge with a specialist may be appropriate Likely to outgrow allergy
α-lactalbumin	+	+	-	
β-lactoglobulin	+	-	-	

*In clinical studies, extensively baked muffin, waffle, and cheese pizza were heated to the point of protein denaturation.

Protein	Bos d 4 / f 76	Bos d 5 / f 77	Bos d 8 / f 78	Management Considerations ^{28,30-38}
α-lactalbumin	+	+	-	<ul style="list-style-type: none"> Avoid all forms of cow's milk Unlikely to become tolerant of cow's milk over time Avoid cow's milk and baked milk products (yogurt, cookies, cakes), as well as products processed with milk (chocolate, sausage, potato chips)
β-lactoglobulin	+	-	-	
Casein	-	+	-	

As in all diagnostic testing, any diagnosis or treatment plan must be made by the clinician based on test results, individual patient history, the clinician's knowledge of the patient, as well as their clinical judgment.

1. Nicolaou N, Poorafshar M, Murray C, et al. Allergy or tolerance in children sensitized to peanut: prevalence and differentiation using component-resolved diagnostics. *J Allergy Clin Immunol.* 2010;125(1):191-197. 2. Bradshaw N, A Clinical Reference Guide to Molecular Allergy. *Go Molecular! Molecular Allergy - The Basics*, 2014. 3. Katelaris CH. Food allergy and oral allergy or pollen-food syndrome. *Curr Opin Allergy Clin Immunol* 2010; 10:246-251. 4. Asarniq A, Nilsson C, Lidholm J, et al. Peanut component Ara h 8 sensitization and tolerance to peanut. *J Allergy Clin Immunol.* 2012;130(2):468-472. 5. Nucera E, et al. Hypersensitivity to major panallergens in a population of 120 patients. *Postepy Dermatol Alergol.* 2015 Aug; 32(4): 255-261. 6. Mittag D, Akkerdaas J, Ballmer-Weber BK, et al. Ara h 8, a Bet v 1-homologous allergen from peanut, is a major allergen in patients with combined birch pollen and peanut allergy. *J Allergy Clin Immunol.* 2004;114(6):1410-1417. 7. Lauer I, Dueringer N, Pokoj S, et al. The non-specific lipid transfer protein, Ara h 9, is an important allergen in peanut. *Clin Exp Allergy.* 2009;39(9):1427-1437. 8. Sastre J. Molecular diagnosis in allergy. *Clin Exp Allergy* 2010; 40:1442-1460. 9. Moverare R, Ahlstedt S, Bengtsson U, et al. Evaluation of IgE antibodies to recombinant peanut allergens in patients with reported reactions to peanut. *Int Arch Allergy Immunol.* 2011;156(3):282-290. 10. Peeters KA, Koppelman SJ, van Hoffen E, et al. Does skin prick test reactivity to purified allergens correlate with clinical severity of peanut allergy? *Clin Exp Allergy.* 2007;37(1):108-115. 11. Asarniq A, Moverare R, Ostblom E, et al. IgE to peanut allergen components: relation to peanut symptoms and pollen sensitization in 8-year-olds. *Allergy.* 2010;65(9):1189-1195. 12. Dang, Thanh, et al. Increasing the accuracy of peanut allergy diagnosis by using Ara h 2. *J Allergy Clin Immunol.* 2012;129(4):1056-1063. 13. Nicolaou N, et al. Quantification of specific IgE to whole peanut extract and peanut components in predication of peanut allergy. *J Allergy Clin Immunol.* 2011;1-2. 14. Vereda A, et al. Peanut allergy. *Clinical and immunologic differences among patients from 3 different geographic regions.* *J Allergy Clin Immunol.* 2010;3(2):1-5. 15. Lange L, Beyer K, Kleine-Tebbe J. Benefits and limitations of molecular diagnostics in peanut allergy. *Allergo J Int.* 2014;23:158-63. 16. Pedrosa, M, et al. Utility of specific IgE to Ara h 6 in peanut allergy diagnosis. *Ann Allergy Asthma Immunol.* 2015;115:108-112. 17. van der Valk, JP, et al. Mono-sensitisation to peanut component Ara h 6: a case series of five children and literature review. *Eur J Pediatr.* 2016;175:1227-1234. 18. Sicherer SH, Sampson HA. Food allergy: Epidemiology, pathogenesis, diagnosis, and treatment. *J Allergy Clin Immunol.* 2014;133(2):291-307. 19. Lemon-Mulé H, Sampson HA, Sicherer SH, Shreffler WG, Noone S, Nowak-Węgrzyn A. Immunologic changes in children with egg allergy ingesting extensively heated egg. *J Allergy Clin Immunol.* 2008;122(5):977-983. 20. Benhamou AH, Caubet JC, Eigenmann PA, et al. State of the art and new horizons in the diagnosis and management of egg allergy. *Allergy.* 2010;65(3): 283-289. 21. Ando H, Moverare R, Kondo Y, et al. Utility of ovomucoid-specific IgE concentrations in predicting symptomatic egg allergy. *J Allergy Clin Immunol.* 2008;122(3):583-588. 22. Shin M, Han Y, Ahn K. The influence of the time and temperature of heat treatment on the allergenicity of egg white proteins. *Allergy Asthma Immunol Res.* 2013;5(2):96-101. 23. Tomicic S, Norman G, Faith-Magnusson K, Jenmalin MC, Devenney I, Böttcher MF. High levels of IgG4 antibodies to foods during infancy are associated with tolerance to corresponding foods later in life. *Pediatr Allergy Immunol.* 2009;20(1):35-41. 24. Urisu A, Yamada K, Tokuda R, et al. Clinical significance of IgE-binding activity to enzymatic digests of ovomucoid in the diagnosis and the prediction of the outgrowing of egg white hypersensitivity. *Int Arch Allergy Immunol.* 1999;120(3):192-198. 25. Bernhisel-Broadbent J, Dintzis HM, Dintzis RZ, Sampson HA. Allergenicity and antigenicity of chicken egg ovomucoid (Gal d III) compared with ovalbumin (Gal d I) in children with egg allergy and in mice. *J Allergy Clin Immunol.* 1994;93(6):1047-1059. 26. Montesinos E, Martorell A, Félix R, Cerdá JC. Egg white specific IgE levels in serum as clinical reactivity predictors in the course of egg allergy follow up. *Pediatr Allergy Immunol.* 2010;21(4 pt 1):634-639. 27. Järvinen KM, Beyer K, Vila L, Bardina L, Mishoe M, Sampson HA. Specificity of IgE antibodies to sequential epitopes of hen's egg ovomucoid as a marker for persistence of egg allergy. *Allergy.* 2007;62(7):758-765. 28. Nowak-Węgrzyn A, Bloom KA, Sicherer SH, et al. Tolerance to extensively heated milk in children with cow's milk allergy. *J Allergy Clin Immunol.* 2008; 122(2): 342-347. 29. Shek LP, Bardina L, Castro R, Sampson HA, Beyer K. Humoral and cellular responses to bovine milk proteins in patients with milk-induced IgE-mediated and non-IgE-mediated disorders. *Allergy.* 2005;60(7):912-919. 30. Wal JM. Bovine milk allergenicity. *Ann Allergy Asthma Immunol.* 2004;93(6 Suppl 3):S2-S11. 31. Sicherer SH, Sampson HA. Cow's milk protein-specific IgE concentrations in two age groups of milk-allergic children and in children achieving clinical tolerance. *Clin Exp Allergy.* 1999;29(4):507-512. 32. Boyano-Martinez T, Garcia-Are C, Pedrosa M, Diaz-Pena JM, Quirce S. Accidental allergic reactions in children allergic to cow's milk proteins. *Allergy Clin Immunol.* 2009;123(4):883-888. 33. Caubet J, et al. Utility of casein-specific IgE levels in predicting reactivity to baked milk. *J Allergy Clin Immunol.* 2013;131(1):222-224. 34. CM Allergy Review. Diagnosis of cow's milk allergy in children: determining the gold standard? *Expert Rev Clin Immunol.* 2014;10(2):257-267. 35. Kim J, et al. Dietary baked milk accelerates the resolution of cow's milk allergy in children. *J Allergy Clin Immunol.* 2011;128(1):125-131. 36. Ito K, et al. The usefulness of casein-specific IgE and IgG4 antibodies in cow's milk allergic children. *Clinical and Molecular Allergy.* 2012;10:1-7.