

ReCAPtulando

Edition n° 57 - Year XIII/2019 - June/July



DUST MITES


Evaluation of sensitization to Der p 1 and Der p 2 in a pediatric population in northern Portugal.

WHEAT

Development of a prediction model for severe wheat allergy.

TRYPTASE

Diagnostic value of Tryptase in Food Allergic reactions: A prospective study of 160 Adult peanut challenges.



Evaluation of sensitization to Der p 1 and Der p 2 in a pediatric population in northern Portugal

Peixoto S, Soares J, Monteiro T, Carvalho M, Santos M, Simões C, et al. *An Pediatr (Brac)* 2018;89(3):162-9.

Commented by:

Prof. Dr. Gustavo Falbo Wandalsen

Master and PhD in Sciences from the Federal University of São Paulo (UNIFESP)

Adjunct Professor of the Discipline of Allergy, Clinical Immunology and Rheumatology of UNIFESP

INTRODUCTION

House dust mites are considered the main source of indoor allergens, and sensitization to them is an important risk factor for the development of respiratory allergic diseases. Among mites, *Dermatophagoides pteronyssinus* (Dp) is the one with the greatest global relevance. About 20 Dp allergenic components have been identified, with Der p 1 and Der p 2 being considered major allergens due to the high sensitization frequency.

Nevertheless, data on the sensitization profile of mite components are poorly understood in many countries. Identification of this profile is of great importance for Dp specific immunotherapy.

METHODS

This is a retrospective study with children and teenagers with asthma and/or allergic rhinitis attended at an allergology service in the northern region of Portugal in 2016. All of them were submitted to an immediate-reading skin test and, in those who tested positive for Dp, ImmunoCAP™ test for Der p 1* and Der p 2² was performed. Specific IgE levels were classified into six ranges: low (0.35 to 0.69 kU_A/L), moderate (0.70 to 3.49 kU_A/L), high (3.50 to 17.49 kU_A/L), very high I (17.50 to <50.00 kU_A/L), very high II (50.00 to <100.00 kU_A/L) and very high III (≥100.00 kU_A/L).

RESULTS

A total of 279 children and teenagers aged 4 to 17 years (65% males) were evaluated, with 71% diagnosed asthma and 88% allergic rhinitis.

Sensitization to Der p 1 was found in 73.5% of patients and Der p 2 in 76.3%. Very high levels (3 categories) to Der p 1 were observed in 41% of cases and Der p 2 in 43%. Only 29 (10%) patients were not sensitized to the two components. Specific IgE levels for Der p 1 and Der p 2 presented a moderate correlation with each other ($r: 0.62$; $p < 0.001$), but poor correlation with Dp papule diameter observed in the skin test (Der p 1 - $r: 0.32$, Der p 2 - $r: 0.25$).

DISCUSSION AND COMMENTS

Mites are the main global source of allergic sensitization and it is estimated that between 65 and 130 million people are sensitized in the world. In this way, a detailed

identification of sensitization profile to these allergens is justified. The present study corroborated the relevance of two main allergenic components of Dp (Der p 1 and Der p 2), which, together, presented allergic sensitization in about 90% of study population.

Changes and home care to reduce indoor exposure and specific immunotherapy are usually indicated for patients with respiratory allergic diseases and mite sensitization. Immunotherapy is the only form of treatment with a potential for change in natural course of allergic diseases, and a reduction in symptoms (nasal and pulmonary) and in the need for rescue medication, as well as improvement of lung function and prevention of new sensitization have been already documented.

Immunotherapy is a long-term, cost-effective treatment and should be prescribed according to individual profile of sensitization and allergens clinical relevance. In Portugal, immunotherapy extracts are standardized only for Der p 1 and Der p 2. Thus, according to the results of this study, about 10% of patients could have mispronounced indication of immunotherapy to Dp, if this was prescribed only by patients results of immediate-reading skin test.

In addition to Der p 1 and Der p 2, other components of Dp are clinically relevant. Der p 10 is a component composed of Dp tropomyosin and cross-reacts with other tropomyosins, such as cockroaches and shrimp. More recently, Der p 23 has been identified as an immunodominant Dp allergen, together with Der p 1 and Der p 2. Studies indicate that Der p 23 plays a relevant role in allergic diseases, especially in asthma.





How to use components in dust mites allergy?

Pedro, 12 years old

BACKGROUND

- Asthma and allergic rhinitis diagnosed several years ago and in regular use of inhaled and nasal corticosteroids.
- Despite environmental measures to reduce exposure to mite allergens, Pedro maintains frequent nasal symptoms and occasional exacerbations of asthma.
- Parents worry about prolonged use of medications.

ALLERGIST EVALUATION

The patient was referred to an allergist who, after analyzing the case, decided to initiate specific immunotherapy.

An immediate response skin test was performed, that was positive for *Dermatophagoides pteronyssinus*, *Blomia tropicalis* and cockroach mix. Besides the skin test, the physician also order the ImmunoCAP[™] components test.

Can allergen components of mites help the proper prescription of immunotherapy?

EXAMS

- ImmunoCAP[™] components test for Der p 1 (d202)* and Der p 2 (d203)*², Der p 10 (d205)*⁵
- ImmunoCAP[™] test for *Blomia tropicalis**⁶

1st scenario

Der p 1 or Der p 2 positive and *Blomia tropicalis* positive: concomitant sensitization to *Dermatophagoides pteronyssinus* and *Blomia tropicalis*, documented by positivity to major components of both mites. Consider concomitant immunotherapy for both mites.

2nd scenario

Der p 1 and Der p 2 negative, Der p 10 positive and *Blomia tropicalis* positive: sensitization only to *Blomia tropicalis*. Possible cross-reaction with other tropomyosins in the skin test for *Dermatophagoides pteronyssinus* (perhaps by sensitization to cockroach). Consider immunotherapy only for *Blomia tropicalis*.

3rd scenario

Der p 1 or Der p 2 positive, *Blomia tropicalis* negative: sensitization to *Dermatophagoides pteronyssinus*, but without evidence of sensitization to *Blomia tropicalis*. Consider immunotherapy only for *Dermatophagoides pteronyssinus*.

Recommended test profile

ImmunoCAP ALLERGENS

Dermatophagoides pteronyssinus (d1)*³
Dermatophagoides farinae (d2)*⁴, *Blomia tropicalis* (d201)*⁶

ImmunoCAP COMPONENTS

Der p 1 (d202)* / Der p 2 (d203)*²

Der p 10 (d205)*⁵
Mite tropomyosin

Clinical Implications

- Allergy-Specific Markers to house dust mites. ⁴
- Patient can benefit from immunotherapy.
- Distinction between sensitization to Der p 1 and Der p 2 help to choose adequate immunotherapy.
- High cross-reactivity between the response to Der p 1 and Der f 1 Der p 2 e Der f 2 do *D. Farinae*. ⁵

- Test helps assess reactivity between the dust mites domestic, crustaceans, insects and molluscs
- Mite tropomyosins such as Der p 10 and Der f 10 have cross reactivity between invertebrates. ⁵
- If the sensitization to Der p 10 is there may be a suspicion of food allergy and patients should be tested against others relevant food components.



Development of a prediction model for severe wheat allergy.

Sugiura S, Matsui T, Furuta T, Sasaki K, Kando N, Ito K. *Pediatric Allergy and Immunology*, v. 29, n. 1, p. 93-96, 2018.

Commented by:

Dra. Ekaterini Goudouris

Professor of the Department of Pediatrics, Faculty of Medicine, UFRJ.

Specialist in Allergy and Immunology by the Brazilian Society of Allergy and Immunopathology

INTRODUCTION

Nowadays the Oral Food Challenge (OFC) has evolved from a definitive diagnostic tool in food allergy to a strategy in determining the food allergen limit dose to initiate oral immunotherapy or even to determine the minimum amount of consumption. However, the OFC test represents a risk of anaphylaxis in patients with severe food allergy.

In order to prevent serious reactions before the OFC, a group of Japanese authors have previously proposed a predictive model for severe egg and cow's milk reactions based on multiple factors.

Wheat is the third antigen most commonly related to food allergy in Japan. Wheat challenge test was identified as an independent risk factor for anaphylaxis at low doses of antigen, requiring the use of adrenaline. Therefore, in a OFC, it is more important to predict severe cases with wheat than it is for other food allergens.

The diagnosis of wheat allergy based on components (omega 5-gliadin) has been shown to be relevant in predicting positive test challenge results but has failed to predict the allergen limit dose or severity of symptoms⁶.

OBJECTIVE

The objective of the present study was to identify the clinical factors that contribute to the severity of the reactions, as well as to develop a predictive model to wheat allergy.

METHODS

OFC tests with udon-type pasta (wheat protein 26mg/g) were performed according to the Japanese guidelines for food allergy of 2014. The severity of the symptoms triggered in the challenge test was determined by the Aichi Anaphylaxis Scoring System (ASCA). They also used the ratio between number of points in ASCA and the cumulative dose of wheat that triggered the symptoms. The diagnosis of anaphylaxis was based on guidelines from the World Allergy Organization (WAO) of 2013.

Development data were obtained from challenge tests with udon-type noodles between April 2012 and December 2014. During this period, 360 tests were performed, 65% of which were positive. Of these, 43 tests were excluded due to lack of complete data.

To validate the results obtained with these data, 123 consecutive provocation tests were prospectively analyzed between January 2015 and March 2016, of which

84 were positive.

The factors that were analyzed as contributors to the OFC test result were: age, gender, history of wheat anaphylaxis, history of atopic dermatitis or asthma, total IgE level, wheat-specific IgE level and omega 5-gliadin-specific IgE level. Patients were classified as severe and non-severe according to the ratio between ASCA points and cumulative protein dose of the allergen that provoked the symptoms, based on the authors' clinical practice.

To compare the sequence data they used a univariate analysis with Mann-Whitney U-test and to analyze the binary variables they used a chi square test. All factors with statistical significance in the univariate analysis were used in the multivariate regression analysis, step by step, to identify independent predictors of severe reactions. Based on the results of these analyzes, the authors constructed a logistic regression model using independent factors with $P < 0.05$. Next, they built a simple model with points. The discriminatory capacity of this model was evaluated using the area under curve ROC. The suitability of the model was tested by the Hosmer-Lemeshow test.

RESULTS

The characteristics of the development group were: average age of 3 years, 68.1% of males. The average ASCA score was 10 with mean cumulative dose of 8 g pasta, resulting in a ratio of 87. Serious cases corresponded to 61% of the cases. Univariate analysis identified 5 significantly different variables among the groups of patients with severe and non-serious reactions: IgE for wheat, IgE for omega 5-gliadin, history of anaphylaxis, age and total IgE level. The multivariate analysis identified as independent factors associated with the severity of the reactions: IgE class for wheat, history of anaphylaxis, total IgE < 1000 kU/L and IgE for omega 5-gliadin positive (> 0.34 kU_A/L). The simple points model constructed has a maximum score of 10 and is shown in table 1 below:

Table 1 - Scoring model to identify risk of serious reactions to wheat

Variables	Number of points
IgE for wheat/class	1/class
IgE for omega 5 gliadin	1
Total IgE < 1000 kU/L	1
Anaphylaxis history	2



How to use components in wheat allergy?

The scoring model created has a significant correlation with the ratio between ASCA points and cumulative dose of wheat that triggers symptoms and was able to predict severity significantly higher than the use of only the wheat-specific IgE class.

DISCUSSION/CONCLUSION

The authors conclude that they developed a simple scoring model to predict serious reactions to the wheat challenge test. The specific IgE for wheat showed a better correlation with the ratio between ASCA points and the protein dose that trigger the symptoms than omega 5-gliadin specific IgE. This may suggest that although omega 5-gliadin is an important component, there are other gliadins, glutenins and water-soluble proteins that may be the cause of allergic reactions.

The fact that a total IgE below 1000 I/mL is an independent risk factor is possibly due to the fact that higher IgE values may promote competition in binding to IgE receptors on the surface of mast cells. This scoring system should not be used in patients with celiac disease because they were excluded from the study.

The proposed scoring system may be clinically useful. Patients with scores ≥ 8 by this system should not be subjected to wheat challenge tests because of the risk of serious reactions. In patients with ≤ 4 , the challenge tests are safe.

COMMENTS

The work brings a simple model, easy to apply and relevant to daily clinical practice. In the daily routine, provocation tests are essential to confirm the diagnosis of food allergy and to plan oral immunotherapy. However, they imply significant serious reactions that we must be prepared and which we should avoid whenever possible. The proposed model includes easily obtainable variables, which include data from the clinical history and complementary tests that are available in our country.

Case Report

Luan, 1 years old

BACKGROUND

- Received exclusive breastfeeding up to the fourth month of life. From that age on, he began to receive formula of

cow's milk, and then there was the introduction of fruits and salt potatoes, without interurrences.

- At 7 months of age, the aunt offered French bread with butter and cheese that was eating in the snack. Approximately one hour later, the patient had very pruritic erythematous plaques on the face, which extended to the trunk and limbs without edema or other clinical changes.

PHYSICIAN EVALUATION

Considering that cow's milk was introduced earlier and without problems, is suspicious of a reaction to wheat. It was decided to carry out a provocation test to clarify the diagnosis, and the suggestion was to carry out the wheat test first. Relatives are quite fearful of taking the test, despite the proposal to perform in a hospital setting.

EXAMS

- ImmunoCAP™ test for Milk (f2)*, casein (f78)*2, beta-lactoglobulin (f77)*3 and alpha-lactalbumin (f76)*4 all below 0,1 kU_A/L
- ImmunoCAP test for wheat (f4)*5 = 2,45 kU_A/L and Tri a 19 (omega 5 gliadin) (f416)*6 = 1,78 kU_A/L

COMMENTS


Patient, according to the suspicion, did not present positive results for cow's milk. The food was already being consumed without problems since the four months of life. The positive result of IgE for wheat confirmed the pediatrician hypothesis of wheat allergy.

We know that the presence of IgE to the omega 5-gliadin component of wheat is related to the risk of positive challenge test. The presence of higher values of specific IgE for wheat, in addition to the clinical history, presence of total IgE <1000 and presence of IgE to omega 5-gliadin should be considered to evaluate the risk of a serious reaction occurring during the test. provocation with wheat.

The patient's tests revealed a positive result for wheat and for wheat omega 5-gliadin, with total IgE <1000 kU/L.

CONCLUSION

In view of the results of the patient's exams, in support of the clinical history, it is possible to say that the patient is allergic to wheat, with a risk of severe reaction during a challenge test, which is therefore not indicated at the present time.



Diagnostic value of Tryptase in Food Allergic reactions: A prospective study of 160 adult peanut challenges.

Dua S, Dowe J, Foley L, Islam S, King Y, Ewan P, Clark AT. *J Allergy Clin Immunol Pract.* 2018 Sep - Oct;6(5): 1692 - 1698. e1. doi: 10.1016/j.jaip.2018.01.006.

Commented by:

Prof. Dr. Marcelo Aun

Medical Alergologist graduated by FMUSP

Professor of the Discipline of Clinical Immunology and Allergy of FMUSP

INTRODUCTION

Serum tryptase is a marker of mast cell activation and an increase in tryptase is a useful indicator of the occurrence of an immediate hypersensitivity reaction. Tryptase has been shown to be valuable in allergic reactions triggered by insect venom and by drugs. However, its usefulness in food allergy remains unknown. Some studies report that tryptase does not increase in food-induced anaphylaxis. Most of these studies rely on data derived from post-mortem samples from patients who died from anaphylaxis or patients who presented acute symptoms in the emergency resulting in varying sampling times and a tendency for more severe reactions. Experimental evidence of tryptase dosage in allergic reactions to food and, in particular, in its typical course of time, are scarce.

OBJECTIVES

The aim of this study was to determine if tryptase increased in allergic reactions to food. In addition, it was evaluated if there would be ideal moments for sample collection and if there would be a diagnostic cutoff point to confirm a clinical reaction.

Results included data from 919 patients, out of 2,369 evaluated, 65.1% male, with a mean age of 3.2 years. A little over 20% had a history of egg-related anaphylaxis, 33.5% had anaphylaxis to other foods. Atopic dermatitis was the allergic disease most frequently associated with egg protein allergy (47.8%), followed by asthma (24%) and allergic rhinitis rhinitis (9.1%). Total IgE levels ranged from 106 to 1040 U/L.

METHODS

Adult patients allergic to peanut were recruited and submitted to up to 4 challenges (oral food challenge [OFC]) with peanut and 1 with placebo, respectively. Tryptase was measured in series on days of challenge before (baseline) and during challenge. The peak percentage of tryptase increasing (peak/baseline) was related to reaction's severity. Receiver Operating Characteristic (ROC) curves were generated by establishing an optimal diagnostic cutoff point.

RESULTS

Tryptase was analyzed in 160 cases with positive OFC (9% with anaphylaxis) and 45 cases with negative

OFC (placebo) tests in 50 adults between 18 and 39 years. Tryptase increased above the normal range (11.4 ng/mL) in only 4 of 160 reactions. However, when compared to basal levels, tryptase increase was observed in 100 of 160 (62.5%) reactions, but in 0 out of 45 placebo challenges. Median of increase (95% confidence interval [CI]) for all reactions was 25% (13.3% to 33.3%) and 70.8% (33.3% to 300%) during anaphylaxis. Peak levels occurred 2 hours after OFC and correlated with severity ($p < 0.05$). Moderate to severe respiratory symptoms, generalized erythema, dizziness and hypotension were correlated with a higher peak/basal tryptase ($p < 0.05$). Analysis of ROC curve showed the optimal cutoff point to identify a reaction as a 30% increase (sensitivity of 0.53, specificity of 0.85), area under the curve of 0.72 (95% IC, 0.67-0.78).


CONCLUSION

Serum tryptase dosage is valuable in allergic reactions to food and correlates with symptoms severity. Comparing peak reaction levels after 2 hours with basal tryptase is essential. An increase in tryptase of 30% is associated with food allergic reactions.

COMMENTS

It's already well defined in literature that there is an association between mastocytosis and anaphylaxis to venoms of hymenoptera and medications. However, this association is not well documented in case of food reactions, nor whether in food anaphylaxis an increase in serum tryptase dosed in acute phase would occur immediately after the reaction.

This study showed that, from the clinical point of view, tryptase levels rarely exceed normal values after an immediate food reaction (4 cases in 160 reactions). However, when compared to basal tryptase, this study showed that tryptase increase occurs after allergic peanut reactions in adults (100 out of 160), and these levels are higher in patients with more severe anaphylactic reactions. In addition, it has been shown that 30% increase in tryptase is associated with food-induced reactions, in this case, peanut. And data showed that 2 hours after exposure to the causative agent was the ideal interval for collection of blood sample.



How to use tryptase in food allergy reactions?

Case Report

Renato, 9 years old

BACKGROUND

• Has moderate to severe intermittent allergic rhinitis since age 6. A history of angioedema at age 8, occurred during a children's buffet party, in which he had eaten several snacks and chocolate candies, with and without nuts. No reaction ever since.

EMERGENCY CARE

Get to the pediatric emergency room with disseminated urticaria, coughing and diarrhea minutes after eating peanuts at his club's bar. Evaluated by the pediatrician, who diagnosed "allergic reaction" and medicated with intravenous corticosteroid and antihistamine. After 15 minutes it evolved with saturation drop (86%) and hypotension. Taken to the emergency room, where received adrenaline 0.3mg intramuscularly in the thigh, with case resolution. An allergist on call at a distance was contacted, which guided collect of serum tryptase, which was performed 1h40min after the start of the reaction. Tryptase value was 9.2 ng/mL (normal to 11.4 ng/mL).

After discharge from the ER, his mother sought the allergist who guided the case and asked for new laboratory tests, collected 50 days after the reaction.

EXAMS

- ImmunoCAP™ Total IgE test*: 0,35kU/L (normal)
- ImmunoCAP Specific IgE for Peanut (f13)*² = 0,35 kU_A/L
- ImmunoCAP Specific IgE for component Ara h 2 (f423)*³ = 14,8 kU_A/L
- ImmunoCAP Tryptase* test: 6,3 ng/mL

Diagnostic value of Tryptase in Food Allergic reactions:

The case illustrates initial reaction situation a year earlier, uninvestigated and that allowed neglect and occurrence of a second, systemic and severe reaction. In the second reaction, the child had an anaphylaxis, unidentified and sub-treated at the beginning, with aggravation in acute phase, but reversed with adequate treatment.

Tryptase dosage close to 2 hours after the reaction did not show a value above the normal, but a new collect after 2 weeks showed that basal tryptase was inferior, that is, the 46% increase in tryptase value confirmed that it was treated with massive degranulation of mast cells, causing anaphylaxis, in this case, a food one. To complete the investigation, the presence of the antigen (peanut and its allergens) has been confirmed and allows the physician to give adequate guidance of future evicton and action plan in case of future inadvertent reactions.

* ImmunoCAP Total IgE 2* ImmunoCAP Allergen f13, Peanut 3* ImmunoCAP Allergen f423, Allergen component rAra h 2 Peanut, 4* ImmunoCAP Tryptase

List of full products names - Main Allergens

ImmunoCAP Allergen f96, Avocado; ImmunoCAP Allergen f210, Pineapple; ImmunoCAP Allergen f225, Pumpkin; ImmunoCAP Allergen f47, Garlic; ImmunoCAP Allergen f255, Plum; ImmunoCAP Allergen f20, Almond; ImmunoCAP Allergen f13, Peanut; ImmunoCAP Allergen f9, Rice; ImmunoCAP Allergen f40, Tuna; ImmunoCAP Allergen f7, Oat; ImmunoCAP Allergen f17, Hazel nut; ImmunoCAP Allergen f92, Banana; ImmunoCAP Allergen f35, Potato; ImmunoCAP Allergen f300, Goat milk; ImmunoCAP Allergen f93, Cacao; ImmunoCAP Allergen f24, Shrimp; ImmunoCAP Allergen f23, Crab; ImmunoCAP Allergen f202, Cashew nut; ImmunoCAP Allergen f18, Brazil nut; ImmunoCAP Allergen f43, Onion; ImmunoCAP Allergen f31, Carrot; ImmunoCAP Allergen f242, Cherry; ImmunoCAP Allergen f6, Barley; ImmunoCAP Allergen f36, Coconut; ImmunoCAP Allergen f340, Codineal extract (Carmine red); ImmunoCAP Allergen f12, Pea; ImmunoCAP Allergen f214, Spinach; ImmunoCAP Allergen f15, White bean; ImmunoCAP Allergen f93, Chicken; ImmunoCAP Allergen f10, Sesame seed; ImmunoCAP Allergen f79, Gluten; ImmunoCAP Allergen f84, Kiwi fruit; ImmunoCAP Allergen f80, Lobster; ImmunoCAP Allergen f33, Orange; ImmunoCAP Allergen f2, Milk; ImmunoCAP Allergen f208, Lemon; ImmunoCAP Allergen f58, Pacific squid; ImmunoCAP Allergen f49, Apple; ImmunoCAP Allergen f293, Papaya; ImmunoCAP Allergen f91, Mango; ImmunoCAP Allergen f247, Honey; ImmunoCAP Allergen f87, Melon; ImmunoCAP Allergen f37, Blue mussel; ImmunoCAP Allergen f8, Maize; ImmunoCAP Allergen f44, Strawberry; ImmunoCAP Allergen f256, Walnut; ImmunoCAP Allergen f245, Egg; ImmunoCAP Allergen f1, Egg white; ImmunoCAP Allergen f75, Egg yolk; ImmunoCAP Allergen f3, Fish (cod); ImmunoCAP Allergen f94, Pear; ImmunoCAP Allergen f95, Peach; ImmunoCAP Allergen f280, Black pepper; ImmunoCAP Allergen f279, Chili pepper; ImmunoCAP Allergen f263, Green pepper (unripe seed); ImmunoCAP Allergen f218, Paprika, Sweet pepper; ImmunoCAP Allergen f26, Pork; ImmunoCAP Allergen f81, Cheese, cheddar type; ImmunoCAP Allergen f82, Cheese, mold type; ImmunoCAP Allergen f14, Soybean; ImmunoCAP Allergen f41, Salmon; ImmunoCAP Allergen f259, Pichard; ImmunoCAP Allergen f259, Grape; ImmunoCAP Allergen f27, Beef; ImmunoCAP Allergen f8, Chicken; ImmunoCAP Allergen f261, Phocodine; ImmunoCAP Allergen f71, Insulin bovine; ImmunoCAP Allergen f70, Insulin porcine; ImmunoCAP Allergen f73, Insulin human; ImmunoCAP Allergen f260, Morphine; ImmunoCAP Allergen f1, Penicillin G; ImmunoCAP Allergen f2, Penicillin V; ImmunoCAP Allergen f85, Alternaria alternata; ImmunoCAP Allergen f33, Aspergillus fumigatus; ImmunoCAP Allergen f5, Candida albicans; ImmunoCAP Allergen f2, Cladosporium herbarum; ImmunoCAP Allergen f209, Penicillium glabrum; ImmunoCAP Allergen f80, Staphylococcal enterotoxin A; ImmunoCAP Allergen f81, Staphylococcal enterotoxin B; ImmunoCAP Allergen f223, Staphylococcal enterotoxin C; ImmunoCAP Allergen f206, Cockroach, American; ImmunoCAP Allergen f6, Cockroach, German; ImmunoCAP Allergen f70, Fire ant; ImmunoCAP Allergen f204, Horse fly; ImmunoCAP Allergen f71, Mosquito; ImmunoCAP Allergen f1, Honey bee venom; ImmunoCAP Allergen f4, Paper wasp venom; ImmunoCAP Allergen f9, Olive; ImmunoCAP Allergen f16, White pine; ImmunoCAP Allergen f18, Eucalyptus, Gum-tree; ImmunoCAP Allergen f119, Acacia; ImmunoCAP Allergen f2, Bermuda grass; ImmunoCAP Allergen f5, Rye-grass; ImmunoCAP Allergen f6, Timothy; ImmunoCAP Allergen f9, Meadow grass, Kentucky blue; ImmunoCAP Allergen f10, Johnson grass; ImmunoCAP Allergen f17, Bahia grass; ImmunoCAP Allergen f201, House dust mite; ImmunoCAP Allergen f2, House dust mite; ImmunoCAP Allergen f3, House dust mite; ImmunoCAP Allergen f352, Allergen component rAra h 8 PR-10; ImmunoCAP Allergen f443, Allergen component rAra o 3, Cashew nut; ImmunoCAP Allergen f354, Allergen component rBer e 1, Brazil nut; ImmunoCAP Allergen f428, Allergen component rCor a 1 PR-10; ImmunoCAP Allergen f439, Allergen component rCyp d 1 Carp; ImmunoCAP Allergen f426, Allergen component rGad cl 1 Cod; ImmunoCAP Allergen f351, Allergen component rHen a 1 Tropomyosin; ImmunoCAP Allergen f419, Allergen component rPru p 1 PR-10; ImmunoCAP Allergen f420, Allergen component rPru p 3 LTP; ImmunoCAP Allergen f421, Allergen component rPru p 4 Pr olin, Peach; ImmunoCAP Allergen f204, Allergen component rBos d 6 BSA, Cow; ImmunoCAP Allergen f222, Allergen component rSus s FSA, Swine; ImmunoCAP Allergen f80, Allergen component rCan f 1 Dog; ImmunoCAP Allergen f102, Allergen component rCan f 2 Dog; ImmunoCAP Allergen f221, Allergen component rCan f 3 Dog serum albumin; ImmunoCAP Allergen f226, Allergen component rCan f 5, Dog; ImmunoCAP Allergen f84, Allergen component rFel d 1 Cat; ImmunoCAP Allergen f220, Allergen component rFel d 2 Cat serum albumin; ImmunoCAP Allergen f228, Allergen component rEqu c 1, Horse; ImmunoCAP Allergen f215, Allergen component rHes v 5 Common wasp; ImmunoCAP Allergen f218, Allergen component rAsp f 1 Aspergillus fumigatus; ImmunoCAP Allergen f229, Allergen component rAlt a 1, Alternaria alternata; ImmunoCAP Allergen f215, Component rGal-alpha-1,3-Gal (alpha-Gal) Throglobulin, bovine; ImmunoCAP Allergen f214, Allergen component MUXF3 CCD, Bromelain

All allergens above refers to ImmunoCAP products. You can find the legal products names on the page before.

EXTRACTS				COMPONENTS		
FOOD	Code	DRUGS	Code	MITES	Code	Functionality
Avocado	f96	Chlorhexidine	c8	Der p 1 (D. pteronyssius)	d202	Cross Reactivity
Pineapple	f210	Folcodine	c261	Der p 2 (D. pteronyssius)	d203	Cross Reactivity
Pumpkin	f225	Bovine Insulin	c71	Der p 10, Tropomyosin (D. pteronyssius)	d205	Cross Reactivity
Garlic	f47	Human insulin	c73	FOOD		
Plum	f255	Swine Insulin	c70	Alfa lactalbumin (Milk)	f76	Species - Specific
Almond	f20	Morphine	c260	Beta-lactoglobulin (Milk)	f77	Species - Specific
Peanut	f13	Penicillin G	c1	Casein (Milk)	f78	Species - Specific
Rice	f9	Penicillin V	c2	Conalbumin (Egg)	f323	Species - Specific
Tuna fish	f40	MICRO-ORGANISMS/FUNGI		Ovalbumin (Egg)	f232	Species - Specific
Oats	f7	Alternaria alternata	m6	Ovomucoid (Egg)	f233	Species - Specific
Hazelnut	f17	Aspergillus fumigatus	m3	Gliadin (Wheat)	f98	Species - Specific
Banana	f92	Aspergillus niger	m207	Omega-5 Gliadin, Tri a 19 (Wheat)	f416	Species - Specific
Potato	f35	Candida albicans	m5	Tri a 14, LTP (Wheat)	f433	Species - Specific
Goat milk	f300	Cladosporium herbarum	m2	Beta-conglycinin (Soy)	f431	Species - Specific
Cocoa	f93	Penicillium notatum	m1	Glicin (Soy)	f432	Species - Specific
Shrimp	f24	Penicillium glabrum	m209	PR-10 (Soy)	f353	Cross Reactivity
Crab	f23	S. enterotoxin A	m80	Ara h 1 (Peanut)	f422	Species - Specific
Cashew nut	f202	S. enterotoxin B	m81	Ara h 2 (Peanut)	f423	Species - Specific
Chestnut	f18	S. enterotoxin C	m223	Ara h 3 (Peanut)	f424	Species - Specific
Onion	f48	S. enterotoxin TSST	m226	Ara h 8, PR-10 (Peanut)	f352	Cross Reactivity
Carrot	f31	INSECTS/VENOMS		Ana o 3 (Cashew nut)	f443	Species - Specific
Cherry	f242	Sewage Cheap	i206	Ber e 1 (Brazil nut)	f354	Species - Specific
Barley	f6	Cockroach	i6	Cor a 1 (Hazelnut)	f428	Species - Specific
Coconut	f36	Foot Washer Ant	i70	Cor a 14 (Hazelnut)	f439	Species - Specific
Carmine Red Dye	f340	Mutuca	i204	Cor a 8, LTP (Hazelnut)	f425	Cross Reactivity
Pea	f12	Stilt	i71	Cor a 9 (Hazelnut)	f440	Species - Specific
Spinach	f214	Venom Bee (Apis mellifera)	i1	Jug r 1 (Walnuts)	f441	Species - Specific
White bean	f15	Venom Vespa / Marimbond	i4	Jug r 3, LTP (Walnut)	f442	Cross Reactivity
Chicken meat	f83	GRAS POLLEN		Parvalbumin from Carpa	f355	Cross Reactivity
Sesame	f10	Olea europaea	t9	Parvalbumin of Cod	f426	Cross Reactivity
Gluten	f79	Salix caprea	t12	Shrimp Tropomyosin	f351	Cross Reactivity
Kiwi	f84	Pinus strobus	t16	Pru p 1, PR-10 (Peach)	f419	Cross Reactivity
Lobster	f80	Eucalyptus spp.	t18	Pru p 3, LTP (Peach)	f420	Cross Reactivity
Orange	f33	Acacia longifolia	t19	Pru p 4, Profilin (Peach)	f421	Cross Reactivity
Milk	f2	Melaleuca leucadendron	t21	ANIMAL		
Lemon	f208	TREE POLLEN		Bovine serum Albumin	e204	Cross Reactivity
Squid	f58	Cynodon dactylon	g2	Pig serum albumin	e222	Cross Reactivity
Apple	f49	Lolium perenne	g5	Can f 1 (Dog)	e101	Species - Specific
Papaya	f293	Pheum pratense	g6	Can f 2 (Dog)	e102	Species - Specific
Mango	f91	Poa pratensis	g8	Can f 3 (Dog serum albumin)	e221	Cross Reactivity
Honey	f247	Sorhum halepense	g10	Can f 5 (Dog)	e226	Species - Specific
Melon	f87	Paspalum notatum	g17	Fel d 1 (Cat)	e94	Cross Reactivity
Blue Mussel	f37	DUST MITES		Fel d 2 (Cat serum albumin)	e220	Cross Reactivity
Corn	f8	Blomia tropicalis	d201	Fel d 4 (Cat)	e228	Species - Specific
Strawberry	f44	D. farinae	d2	Equ c1 (Horse)	e227	Species - Specific
Walnuts	f256	D. microceras	d3	LATEX		
Egg	f245	D. pteronyssinus	d1	Hev b 1 (Latex)	k215	Species - Specific
Egg, clear	f1	House dust	h2	Hev b 3 (Latex)	k217	Species - Specific
Egg yolk	f75	ANIMAL		Hev b 5 (Latex)	k218	Cross Reactivity
Fish	f3	Dog	e5	Hev b 6.02 (Latex)	k220	Cross Reactivity
Wait	f94	Horse	e3	Hev b 8, Profilin (Latex)	k221	Cross Reactivity
Peach	f95	Chicken	e85	Hev b 11 (Latex)	k224	Cross Reactivity
black pepper	f220	Cat	e1	VENOMS		
Red pepper	f279	Caw	e4	Api m 1 - Fosfolipase A2 (Bee)	i208	Species - Specific
green pepper	f263	OTHERS		Pol d 5 (Marimbond)	i210	Species - Specific
Pepper	f218	Cotton	o1	Ves v 1 - Fosfolipase A1 (Wasp)	i211	Species - Specific
Octopus	f59	Latex	k82	Ves v 5 (Wasp)	i209	Species - Specific
Pork meat	f26	Tobacco Sheet	o201	FUNGI		
Cheese (type Cam, Brie, Roqf)	f82	MISCELLANEOUS		Asp f 1 (Aspergillus fumigatus)	m 218	Species - Specific
Cheese (Cheddar type)	f81			Alt a 1 (Alternaria alternata)	m 229	Species - Specific
Soy	f14					
Salmon	f41					
Sardine	f61					
Grape	f259					
Wheat	f4					
Black wheat	f11					
Cow meat	f27					

References: 1. Evaluation of sensitization to Der p 1 and Der p 2 in a pediatric population in northern Portugal. Peixoto S, Soares J, Monteiro T, Carvalho M, Santos M, Simões C, et al. An Pediatr (Brac) 2018;89(3):162-9. 2. Development of a prediction model for severe wheat allergy. Sugiura S, Matsui T, Furuta T, Sasaki K, Kando N, Ito K. Pediatric Allergy and Immunology, v. 29, n. 1, p. 93-96, 2018. 3. Diagnostic value of Trypstatin in Food Allergic reactions: A prospective study of 160 adult peanut challenges. Dua S, Dowey J, Foley L, Islam S, King Y, Ewan P, Clark AT. J Allergy Clin Immunol Pract. 2018 Sep - Oct;6(5): 1692 - 1698. e1. doi: 10.1016/j.jaip.2018.01.006. 4. Thomas WR. Human Vaccines & Immunotherapeutics 2012;8(10):1469-78. 5. Bessot JC and Pauli G. Eur Ann Allergy Clin Immunol. 2011;43:141-56. 6. Ito K, Futamura M, Borres MP, et al. IgE antibodies to omega-5 gliadin associate with immediate symptoms on oral wheat challenge in Japanese children. Allergy. 2008;63:1536-1542.