

MMDx Kidney, Heart & Lung

Personalized Transplant Care Through Precision Medicine

Molecular Microscope®
Diagnostic System **MMDx**

What is the Molecular Microscope® Diagnostic System?

MMDx® Kidney, MMDx® Heart and MMDx® Lung are biopsy-based Laboratory Developed Tests that measure gene expression profiling and provide risk assessment for rejection and injury in transplant organs. These tests, which can be used in conjunction with conventional histopathology, measure mRNA transcript levels in the biopsy sample and apply an algorithm to score¹ the results.

The MMDx Kidney, Heart or Lung tests can help stratify the risk for conditions like T-cell mediated rejection (TCMR), antibody-mediated rejection (ABMR), acute and chronic injury, atrophy fibrosis, and arterial hyalinosis.

A New Era of Precision Medicine

MMDx combines gene expression technology with the power of big data to deliver objective and reproducible transplant biopsy assessments.


The MMDx algorithm compares the molecular features of the biopsy to those in an organ-specific reference set. The organ-specific reference sets were developed over the course of a decade and consist of over 1600 kidney samples, 1000 heart samples, 600 transbronchial biopsies (TBB) lung samples, and 300 third airway bifurcation (3BMB) mucosal samples.

The reference sets incorporate data that is representative of the global transplant population from early post-transplant to more than 30 years post-transplant.

Concordance to Histological Biopsy Diagnosis

Histopathology scores are the primary tool for diagnosing injury or rejection, but studies show frequent disagreement in histological TCMR diagnosis.

Additionally, the quality of the biopsy samples can sometimes impact histology results, rendering the samples unreadable. The MMDx tests for Kidney, Heart and Lung are not intended to replace histology. Rather, they can be used in addition to a histopathologist's assessment, especially for the objective evaluation of challenging cases.



Variability in TCMR Diagnosis with Histopathology

When assessing the same heart biopsy sample, only 28% of pathologists agreed on the diagnosis of TCMR². Likewise, a mere 18% agreed on a TCMR diagnosis for a single lung sample³ and 50% were in agreement after reviewing the same kidney sample¹.

Advantages of MMDx vs. Histology

The **MMDx Kidney, Heart and Lung tests complement conventional biopsy** and improve the assessment of **rejection** and **injury** in transplanted organs.

Actionable data: Provides objective and probabilistic risk assessment

Fast turnaround: Results available up to 1-2 business days after receipt of the sample

Easily Incorporated: Sample is taken from the existing biopsy (requires only 3 mm for kidney, 1-2 bites for heart and lung)

Efficient Process: Place biopsied tissue into a tube (prepared with RNAlater®) and ship at room temperature

What Does MMDx Assess?

Rejection related:

- Risk of T-cell mediated rejection (TCMR)
- Risk of Antibody-mediated rejection (ABMR)

Injury & non-rejection related disease:

- Acute and chronic injury
- Atrophy fibrosis
- Arterial hyalinosis

Reducing the Economic Impact of Biopsy Testing

According to a study on the economic impact of molecular and histological biopsy testing⁴, the application of MMDx Kidney along with histological assessment can improve the diagnosis of graft dysfunction and may even produce cost savings in rejection-related treatment.

By more accurately characterizing rejection, the gene expression profiling test may support reductions in costs associated with graft failure, such as hospitalization, dialysis and repeat transplantation.

The study suggests that the clinical use of MMDx Kidney has generated an estimated undiscounted savings of **\$19,271 per test** over a five-year period, assuming one test per patient.

Current and Emerging Applications for Molecular Assessment

Research suggests that molecular assessment may be useful, or even necessary, to accurately estimate rejection and injury in transplanted organs⁵.

MMDx Kidney

In kidney patients with chronic active antibody-mediated rejection (caAMR) and a high degree of chronicity, molecular evidence of rejection has been used to track responses to immunosuppressive therapies and identify response to treatment, as evidenced by improved inflammation⁶.

Patient Information

Clinical Interpretation

Summary of Molecular Changes (injury, rejection, scores)

Visualization Relationship of biopsy to others in reference set

Indicates this biopsy

Molecular Microscope[®] Diagnostic System MMDx

Molecular Microscope[®] Diagnostic Report for Kidney (MMDx-Kidney)

Patient and Institution Information

Patient Name or ID	PC191210 191203009MM	Lab ID Number	191210004MM
Patient DOB	Not Provided	Ordering Physician	Dr. S Jones
Patient MRN	Not Provided	Ordering Institution	Kashi Clinic TEST ACCT

Testing and Clinical Information

Test Date	12/10/2019	Time from Transplant to Biopsy	1 years
Report Date	12/11/2019	DSA Status	Not Provided
Transplant Date	05/26/2018	Biopsy Indication	Surveillance, compare to previous
Biopsy Date	12/02/2019	Primary Disease	HTN

Pure Molecular Interpretation (Results Summary)


Abnormal kidney transplant biopsy. Moderate early-stage ABMR. No TCMR. Mild atrophy-fibrosis with minimal inflammation and AKI. Compared to the biopsy of June 26th, 2019, there has been resolution of inflammation and AKI features. Note: the Molecular Microscope [®] Diagnostic System cannot exclude primary glomerular diseases.	Percent cortex ¹
	96%

Result Details

Biopsy Rejection and Injury Scores

	Classifier / Gene Sets	Biopsy Score	Range of Values ²	Upper Limit of Normal ³	Interpretation
Injury Scores	Global Disturbance Score	-1.19	-3.8 — 5.8	0.03	Minimal
	Acute Kidney Injury (AKI) Score	-0.13	-0.6 — 1.6	0.55	Mild
	Atrophy-Fibrosis Score	0.34	0 — 1	0.52	Mild
	Rejection Score	0.57	0 — 1	0.30	Moderate
Rejection Scores	T Cell-Mediated Rejection (TCMR) Score	0.01	0 — 1	0.10	Normal
	Antibody-Mediated Rejection (ABMR) Score	0.56	0 — 1	0.20	Moderate

Archetypal Analysis (please see Archetypal Analysis Description on Page 2 for details)



Rejection phenotype ⁴ (six scores, R1-R6, adding up to 1.0)			
R1 Non-rejecting	0.13	R4 Early-Stage ABMR (EABMR)	0.62
R2 TCMR	0.00	R5 Fully-Developed ABMR (FABMR)	0.25
R3 Mixed Rejection	0.00	R6 Late-Stage ABMR (LABMR)	0.00
All ABMR (Sum of R4, R5, and R6)			0.87

1. Percent cortex is a quality control measure.

2. The 2.5th to 97.5th percentiles in the entire Reference Set.

3. 90th percentile in relevant Reference Set biopsies.

4. Scores from archetypal analysis.

MMDx Kidney: Sample Report


MMDx Kidney (cont.)

Inflammation in areas of atrophy-fibrosis (i-IFTA) has shown to be associated with increased risk of failure in kidney biopsies. A recent study has concluded that i-IFTA in indication biopsies reflect current or ongoing parenchymal injury, either with TCMR or (more commonly) with concomitant ABMR⁷.

MMDx Kidney has demonstrated accuracy and reproducibility in kidney biopsy assessment with minimal inter-observer disagreement in reporting. As a result, MMDx may be particularly valuable in cases when pathology results are “borderline” or “suspicious”⁸.

MMDx Heart


MMDx has demonstrated an improved sensitivity for the detection of subclinical graft injury in heart patients⁷ and for the detection of ABMR⁷. Molecular assessment may also provide greater clarity for biopsies with high probability of molecular injury but no molecular rejection, a state that can often be misdiagnosed as rejection by histology⁹.



MMDx Heart: Sample Report

MMDx Lung

For lung patients, the histological assessment of transbronchial biopsies (TBB) offers limited reproducibility and presents considerable risk. However, evidence indicates molecular assessment with MMDx is a promising diagnostic solution for mucosal biopsies, a method that is much safer but not histologically interpretable¹⁰.



The MMDx Process



Ordering Information

- Click "Specimen Kit Ordering"
- Complete ordering form and click "Order"
- A Specimen Kit will be delivered to your center



Specimen Kit Includes:

- Instructions
- Specimen tube and labels
- Return containers and labels



Specimen Handling

- Unpack your Specimen Kit and instructions.
- Place biopsy (Kidney: 3-5 mm of biopsy or portion of core) immediately from the needle or bite (Heart and Lung: 2 bites) in specimen tube provided in the kit. Washes will destroy the sample; please refrain from applying formalin, saline, or any other washes.
- Pack and ship. Pack samples and ship to address on the label provided in the kit.



Biopsy Results

- Results are available within 1-2 days after receipt of sample.
- Upon receipt of your sample, an account will be created for you on a secure client portal.
- You will receive an email notification once the report is available on the portal.



Client Services Assistance

- Assist with result interpretation/consultation
- Specimen requirements and handling
- Requisition completion
- Report status
- Client supplies
- Cost and billing

Service Lab & Billing Information

Kashi Clinical Laboratories, Inc.
Portland, Oregon
Phone: (877) 879-1815 or (503) 206-4989
Fax: (503) 206-6939

Hours of Operation

Monday – Friday: 7:00 am – 7:00 pm (PT)
Weekend and Holidays: on-call coverage



Sales & Technical Support

Email: 1lambda-mmdx@thermofisher.com
Website: go.1lambda.com/mmdx

One Lambda | A Thermo Fisher Scientific Brand
22801 Roscoe Blvd, West Hills, CA 91304, USA
TEL: 747.494.1000
FAX: 747.494.1001

International: Contact your local distributor

References

1. Reeve, J. et al. Generating automated kidney transplant biopsy reports using ensembles of molecular classifiers. *Am J Transplant* 19 (10): 2719-2731, 2019
 2. Crespo-Leiro, M. G. et al. Concordance among pathologists in the second Cardiac Allograft Rejection Gene Expression Observational Study (CARGO II). *Transplantation*. 2012 Dec 15;94(11):1172-7. doi: 10.1097/TP.0b013e31826e19e2. PMID: 23222738.
 3. Arcasoy, S. M. et al. Pathologic interpretation of transbronchial biopsy for acute rejection of lung allograft is highly variable. *Am J Transplant*. 2011 Feb;11(2):320-8. doi: 10.1111/j.1600-6143.2010.03382.x. Epub 2011 Jan 10. PMID: 21219569.
 4. Fustfeld, L., Menon, S., Gupta, G., Lawrence, C., Masud, S. F., & Goss, T. F. (2022). U.S. payer budget impact of a microarray assay with machine learning to evaluate kidney transplant rejection in for-cause biopsies. *Journal of Medical Economics*. DOI: 10.1080/13696998.2022.2059221
 5. Halloran, P. F. et al. Exploring the cardiac response to injury in heart transplant biopsies. *JCI Insight*. 2018 Oct 18;3(20):e123674. doi:10.1172/jci.insight.123674. PMID:30333303; PMCID: PMC6237487.
 6. Kumar, D. et al. Impact of Belatacept Conversion on Renal Function, Histology, and Gene Expression in Kidney Transplant Patients With Chronic Active Antibody-mediated Rejection. *Transplantation*. 2021 Mar 1;105(3):660-667. doi: 10.1097/TP.0000000000003278. PMID: 32510913.
 7. Halloran, P. F. et al. Molecular phenotype of kidney transplant indication biopsies with inflammation in scarred areas. *Am J Transplant*. 2019 May;19(5):1356-1370. doi: 10.1111/ajt.15178. Epub 2018 Dec 13. PMID: 30417539.
 8. Madill-Thomsen, K. S. et al. Discrepancy Analysis Comparing Molecular and Histology Diagnoses in Kidney Transplant Biopsies. *Am J Transplant* https://doi.org/10.1111/ajt.15752. *Am J Transplant*. 2020;20: 1314-1350
 9. Nguyen, V. P., & Kobashigawa, J. A.. Antibody-mediated rejection after heart transplantation: diagnosis and clinical implications. *Curr Opin Organ Transplant*. 2020 Jun;25(3):248-254. doi: 10.1097/MOT.0000000000000754. PMID: 32304428.
 10. Halloran, K. et al. Molecular phenotyping of rejection-related changes in mucosal biopsies from lung transplants. *Am J Transplant*. 2020 Apr;20(4):954-966. doi:10.1111/ajt.15685. Epub 2019 Dec 1. PMID: 31679176.
- * Reeve J, et al. *Am J Transplant* 2013.
† Crespo-Leiro MG, et al. *Transplant* 2012.
‡ Arcasoy SM, et al. *Am J Transplant* 2011.

■ To learn how you can incorporate MMDx into your clinical practice, contact us at **1lambda-mmdx@thermofisher.com** or visit **go.1lambda.com/mmdx**.

One Lambda | A Thermo Fisher Scientific Brand
22801 Roscoe Blvd, West Hills, CA 91304, USA
TEL: 747.494.1000
FAX: 747.494.1001
International: Contact your local distributor

MMDx and Molecular Microscope are registered trademarks of Transcriptome Sciences Inc. MMDx Heart, MMDx Kidney, and MMDx Lung are tests that were developed and validated by Kashi Clinical Laboratories, Inc. The laboratory developed tests are used for clinical purposes by the CLIA-certified laboratory performing the test. These tests have not been cleared or approved by the US FDA or CE marked in the EU as an in vitro diagnostic test. RNeasy is a trademark of Qiagen, Inc.
© 2022 Thermo Fisher Scientific. All rights reserved.



A Thermo Fisher Scientific Brand