

L. Y. Wong, P. Vatta, P. Balachandran, R. Tebbs, C. Cummings, P. Brzoska, Manohar Furtado, and Olga Petruskenskaya  
Applied Biosystems, Foster City, CA

### Abstract


*Escherichia coli* O157:H7 poses a major threat to public safety, as evidenced by recent outbreaks of illnesses stemming from contaminated food supplies including ground beef and spinach. *E. coli* O157:H7 causes hemorrhagic colitis and can lead to severe complications that are potentially fatal. There is an apparent need for a rapid, reliable, and specific assay to detect and identify the O157:H7 serotype. The aim of the study was to develop an assay to detect O157:H7 while avoiding detection of other serotypes of *E. coli*, including O55:H7 which is genetically very closely related to O157:H7. Using publicly available sequences of O157:H7 and other *E. coli* serotypes, putative O157:H7 specific sequences were identified. However, genomic information for O55:H7 was extremely limited. To distinguish between O55:H7 and O157:H7, the genome sequence of O55:H7 was determined using the SOLiD™ system sequencing approach and compared to O157:H7. Whole genome comparison of O157:H7 and O55:H7 sequences yielded genome segments unique to O157:H7. These O157:H7 specific sequences were used as potential targets to design TaqMan® real-time PCR assays. The computational assay design pipeline minimized mismatches to O157:H7 sequences while maximizing mismatches to all other exclusion targets. The TaqMan assays were determined to be specific for the O157:H7 serotype by testing a large panel of *E. coli* of various serotypes and related pathogens such as *Shigella*. The assays exhibited no cross-reactivity against overnight cultures of different food matrices which confirmed that the assays did not detect other microorganisms that might be present in the background flora of the samples. High sensitivity of the assays in conjunction with novel sample prep solutions allowed detection of 1 to 3 cfu of *E. coli* O157:H7 in 65 g of cultured ground beef during an 8 hour sample-to-result workflow period. In summary, this study demonstrated that *E. coli* O157:H7 possesses unique sequences that can be targeted for sensitive detection and accurate differentiation of O157:H7 from other closely related *E. coli* serotypes including O55:H7.

### *E. coli* O157:H7-specific sequences

Based on genomic comparisons, we designed several TaqMan real-time PCR assays against the serotype *E. coli* O157:H7 and avoided detection of other serotypes. *E. coli* O55:H7, the nearest neighbor of O157:H7, is extremely close in sequence to O157:H7. Because genomic sequence of O55:H7 was not readily available for genomic comparison prior to our design process, we generated its genomic sequence using Applied Biosystems' SOLiD technology to compare with an O157:H7 reference sequence. Genomic regions that were specific to O157:H7—either different or absent in O55:H7—were identified as putative target sequences for the O157:H7 TaqMan primers and probes.

### TaqMan® and the assay design process

Figure 3. TaqMan assay development pipeline for *E. coli* O157:H7 targets. A TaqMan assay consists of 2 primers and a dye-labeled probe and non-fluorescent quencher.



### Specificity

#	Strain ID	Bacterial strain	Assay 1, Ct	Assay 2, Ct
1	PE30	<i>E. coli</i> O157:H7	33.2	35.3
2	PE677	<i>E. coli</i> O157:H7	29.0	29.5
3	PE678	<i>E. coli</i> O157:H7	28.9	36.1
4	PE679	<i>E. coli</i> O157:H7	31.1	33.6
5	PE700	<i>E. coli</i> O157:H7	33.1	34.9
6	PE701	<i>E. coli</i> O157:H7	32.0	33.0
7	PE702	<i>E. coli</i> O157:H7	34.0	35.0
8	PE703	<i>E. coli</i> O157:H7	34.0	35.5
9	PE40	<i>E. coli</i> O157:H7	28.0	29.4
10	PE474	<i>E. coli</i> O55:H7	No signal	No signal
11	PE704	<i>E. coli</i> O55:H7	No signal	No signal
12	PE705	<i>E. coli</i> O55:H7	No signal	No signal
13	PE706	<i>E. coli</i> O55:H7	No signal	No signal
14	PE735	<i>E. coli</i> O55:H7	No signal	No signal
15	PE684	<i>E. coli</i> O137:H4†	No signal	No signal
16	PE37	<i>E. coli</i> O145:NM	No signal	No signal
17	PE526	<i>E. coli</i> O154:H25	No signal	No signal
18	PE827	<i>E. coli</i> O156:H8	No signal	No signal
19	PE731	<i>E. coli</i> O26:H32	No signal	No signal
20	PE688	<i>E. coli</i> O28:H35	No signal	No signal
21	PE687	<i>E. coli</i> O48:H21	No signal	No signal
22	PE673	<i>E. coli</i> O5:NM	No signal	No signal
23	PE39	<i>E. coli</i> O78:K80:H12	No signal	No signal
24	PE1199	<i>Shigella flexneri</i>	No signal	No signal
25	PE1200	<i>Shigella boydii</i>	No signal	No signal
26	PE1201	<i>Shigella sonnei</i>	No signal	No signal





Figure 4. TaqMan results when 2 different *E. coli* O157:H7 assays were tested on a panel of *E. coli* O157:H7 and non-O157:H7 strains. Ct (cycle threshold) values are indicated. Non-O157:H7 strains including O55:H7 did not show any cross-reactivity with the assays. The reactions were run on the AB 7500 Fast Sequence Detection System.

### Sensitivity




Figure 5. Fig 5a. Assay efficiencies close to 100% were determined for 2 assays by testing the assays over a dilution series of *E. coli* O157:H7 purified genomic DNA. 10 copies were consistently detected by both assays. The internal positive control appears in red. Fig 5b. 1 to 3 cfu of *E. coli* O157:H7 was spiked in 65 g ground beef, enriched for 6 hr, and detected by the 2 assays post-sample prep.

### Sample-to-result workflow





Figure 6. Sample prep is a crucial step in detection of *E. coli* O157:H7 within various food matrices. Foods are complex matrices and, without the appropriate sample prep, can often inhibit PCR. Avoiding inhibition is especially critical in order to enable detection of low levels of bacteria in a food sample. Fig. 6a outlines the steps involved in testing various food matrices for pathogens. Using Applied Biosystems' sample prep and Real-Time detection system, *E. coli* O157:H7 can be enriched in a food sample, prepared, and run on TaqMan with results within an 8-hour workflow. Fig. 6b, depicts one of Applied Biosystems' novel sample prep options. Bacteria in a food matrix is lysed, its nucleic acids separated using magnetic particles, and then eluted in a buffer compatible with PCR.



### Spiked and unspiked food matrices






Figure 7. Spinach, whole milk, and mayonnaise were spiked with 1-3 cfu *E. coli* O157:H7. The cultures were enriched overnight in this study. Following sample prep, two *E. coli* O157:H7 assays were run on each sample (signals shown in blue and purple). Unspiked controls were included. Absence of signal from unspiked samples indicates no cross-reactivity of the assays to the food matrices. The internal positive control (yellow) is included in the matrix. Overnight enrichment of 1-3 cfu *E. coli* O157:H7 was detected by both assays in almost all cases. Only in figure 7c, spike #4 sample did not enrich well for *E. coli* O157:H7 and generated high Ct for both assays.

### Conclusions

- Highly specific TaqMan assays for the detection of *E. coli* O157:H7 were designed using the sequenced genomes of O157:H7 and O55:H7 and applied Applied Biosystems' assay design pipeline. The assays were able to distinguish between *E. coli* O157:H7 and its nearest neighbor O55:H7 in experimental tests.
- The sensitivity of the assays permits reliable detection of 10 copies of *E. coli* O157:H7 genomic DNA. It also detects 1-3 cfu of *E. coli* O157:H7 after enrichment in food matrices such as ground beef, whole milk, and mayonnaise.
- AB sample prep removes PCR inhibitors and permits detection of low level of *E. coli* O157:H7 enriched in complex food matrices.
- The assays do not cross-react with ground beef, spinach, whole milk, or mayonnaise.
- We have developed an 8 hr sample-to-result workflow for ground beef.

### Trademark/Licensing

For Research Use Only. Not for use in diagnostic procedures.

Practice of the patented 5' Nuclease Process requires a license from Applied Biosystems. The purchase of the TaqMan assay includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research when used with the separate purchase of an Authorized 5' Nuclease Core Kit. No other patent rights are conveyed expressly, by implication, or by estoppel. For further information on purchasing licenses contact the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

Applied Biosystems, and AB (Design) are registered trademarks of Applied Corporation or its subsidiaries in the US and/or certain other countries. TaqMan is a registered trademark of Roche Molecular Systems, Inc. All other trademarks are the sole property of their respective owners.