



# Rapid Detection of *Salmonella* in Pet Food: Design and Evaluation of Reduced Time-to-Result workflow Using Real-Time PCR



P.Balachandran<sup>1\*</sup>, M. Brevnov<sup>1</sup>, M. Furtado<sup>1</sup>, O. Petrauskene<sup>1</sup>, R. Tebbs<sup>1</sup>, M. Friberg<sup>2</sup>, V. Vanlandingham<sup>2</sup>, K. Kozak<sup>2</sup> and D. Charbonneau<sup>2</sup>

<sup>1</sup> Life Technologies, 850 Lincoln Centre Drive, Foster City, California

<sup>2</sup> Procter & Gamble Cincinnati Ohio

Contact  
Priya Balachandran, PhD  
Priya.Balachandran@lifetech.com

## INTRODUCTION

According to the 2009/2010 National Pet Owners Survey, 62% of U.S. households own a pet, and the estimated sales within the U.S. market for pet food alone is \$17.4 billion. Consequently, rapid and robust testing of pet food for food borne pathogens is gaining importance. Food samples create unique challenges, particularly with respect to sample preparation for detection by PCR. This is especially true for high carbohydrate containing matrices, such as pet food. To address these needs, we have designed and evaluated a complete and streamlined workflow to detect *Salmonella* in pet food that includes:

- **Short Pre-enrichment** : Pre-enrichment is conducted under shaking conditions to enhance rapid growth of *Salmonella*.
- **Novel, custom-designed sample preparation**: A novel pre-clarification tray was designed for addressing difficult sample matrices like pet food, which when combined with automated sample preparation, enables high-throughput processing.
- **Lyophilized Real-Time PCR assay**: Highly sensitive and specific real-time PCR assay allows for addition of more sample volume, which creates the possibility of testing pooled samples.
- **Custom Designed Software** : Simplifies data analyses by providing plus/minus results, and differentiating viable *Salmonella* from background.
- **Shortened Time-To-Result**: The streamlined sample-to-answer workflow offers a great advantage over existing culture confirmation methods by reducing the time-to-result to under 2h following pre-enrichment. Traditional culture confirmation typically takes an additional 2-4 days.

## MATRICES TESTED


For method comparison, a total of 45 samples for each of the four *Salmonella* strains were analyzed. Each sample set included 5 uninoculated samples, 20 low-level and 20 high-level samples. Target levels of each strain were as follows: 0.2-2 colony forming units (cfu) /25 g for the low level inoculation, 2-10 cfu/25g for the high level inoculation and 0 cfu/25 g for the uninoculated control samples. Each pet kibble was inoculated with a different strain of *Salmonella*. A 24 hour broth culture inoculum was utilized to inoculate a bulk sample of ground pet food and kept for 2 weeks at room temperature (25°C) as per AOAC guidelines. For enumeration a 3-tube MPN analysis of each inoculated food was used to determine the two target inoculum levels. On the day of initiation of analysis, samples of 100 g, 10 g, 1 g and 0.1 g from each inoculum level were prepared for each matrix, enriched with the appropriate medium and confirmed according to the FDA BAM reference method procedures.

Sample Identification	Inoculating Organism	Sample Identification	Inoculating Organism
1-45	<i>S. enterica</i> ser. Typhimurium ATCC 14028	101-145	<i>S. enterica</i> ser. Senftenberg ATCC 43845
201-245	<i>Salmonella enterica</i> subspecies Tennessee	301-345	<i>Salmonella enterica</i> subspecies Othmarschen


For research use only. Not intended for human or animal therapeutic or diagnostic use. The PrepSEQ™ Nucleic Acid Isolation Kit is manufactured and sold under license from GE Healthcare under U.S. Patent Nos. 5,523,231 and 5,681,346 and other foreign patents. End Users are specifically not authorized to and are forbidden from reselling, transferring or distributing any products either as a stand alone product or as a component of another product. © 2010 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. TaqMan is a registered trademark of Roche Molecular Systems, Inc.

## INTEGRATED SYSTEM FOR STREAMLINED SAMPLE-TO-RESULT WORKFLOW

### Pre-Enrichment



### Sample Preparation



### Real-time PCR and Data Analyses



**Acknowledgements:** We would like to thank Karl Meyer, Amanda Manolis and Brooke Schwartz (LIFE Technologies) for their significant contributions.

## Integrated System Performance

Figure 1: Preclarification Improves Signal Intensity by Reducing PCR Inhibitor Load





Figure 2: Validation of Complete Workflow Demonstrated Equivalency between the MicroSEQ® *Salmonella* spp. Method and Reference Method



MicroSEQ® <i>Salmonella</i> spp. vs. FDA-BAM					
Inoculation Level	MPN/25g	MicroSEQ Presumptive +	MicroSEQ Confirmed +	FDA-BAM Confirmed +	χ <sup>2</sup> (MicroSEQ vs. FDA-BAM)
<b><i>S. enterica</i> ser. Typhimurium ATCC 14028</b>					
Control	0	0/5	0/5	0/5	-
Low	0.1	0/20	0/20	0/20	-
High	2.3	6/20	5/20	6/20	0.12
<b><i>S. enterica</i> ser. Senftenberg ATCC 43845</b>					
Control	0	0/5	0/5	0/5	-
Low	0.1	4/20	3/20	3/20	-
High	0.4	12/20	13/20	14/20	0.43
<b><i>S. enterica</i> ser. Tennessee</b>					
Control	0	0/5	0/5	0/5	-
Low	2.3	2/20	2/20	6/20	2.44
High	11.5	18/20	18/20	19/20	0.35
<b><i>S. enterica</i> ser. Othmarschen</b>					
Control	0	0/5	0/5	0/5	-
Low	0.2	6/20	6/20	7/20	0.11
High	2.3	16/20	16/20	18/20	0.76

\*Values > 3.84 are significant at 5%

## Conclusions

- Inclusion of pre-clarification during sample preparation addressed problems due to matrix carryover, including PCR inhibition.
- The entire workflow can be completed in a shortened time frame.
- Inclusion of software analyses enabled differentiation of live *Salmonella* from background.
- Validation studies demonstrated equivalency between the presented method and the reference FDA-BAM protocol.