

Highly Specific and Sensitive Detection and Quantification of *Staphylococcus aureus* Using Real-Time PCR



Lily Wong¹, Pius Brzoska¹, Timothy Dambaugh², Dawn Fallon², Morgan Wallace², Bridget Andaloro², George Tice², Manohar Furtado¹, and Olga Petruskenskaya¹

¹ Applied Biosystems, Foster City, CA, ² DuPont Qualicon, Wilmington, DE

Abstract

Staphylococcus aureus is a Gram-positive pathogen that can contaminate a variety of foods. *S. aureus* is found in meat, poultry, eggs, milk and dairy products, salads such as egg and tuna, and baked goods with cream fillings. Specific and sensitive detection of *S. aureus* has been achieved using an assay based on careful target selection and bioinformatics algorithms.

Purpose: The aim of this study was to develop a specific, reliable, and simple-to-use Real-Time PCR assay for *S. aureus* and demonstrate its utility for rapid and sensitive screening of enrichments of meat and dairy products.

Methods: A duplex Real-Time PCR assay was developed to detect *S. aureus* and a positive control sequence to monitor PCR inhibition. Assay performance was tested on an inclusivity panel of *S. aureus* isolates and an exclusivity panel of related *Staphylococcal* species. Sensitivity of the assay was evaluated using a quantified *S. aureus* DNA sample.

Results: The *S. aureus* assay, tested in lyophilized and tableted form, detected 27 out of 27 (100%) *S. aureus* isolates in the inclusivity panel and none of the 25 (0%) species in the exclusivity panel. Limit of detection was estimated at 10 genomic copies of *S. aureus*. In naturally contaminated ground beef, *S. aureus* levels of 10 cfu per gram were detected within 24 hrs.

Significance: This study demonstrated that our assay specifically detects *S. aureus* and enables accurate identification of *S. aureus* from other *Staphylococcal* species. This specificity combined with the high sensitivity of the assay and ease-of-use would be amenable to *S. aureus* screening of diverse foods.

Assay design and features

The primers and probes in the *S. aureus* Real-Time PCR (TaqMan®) assay were designed using Applied Biosystems' rigorous bioinformatics pipeline that tabulates the effects of potential mismatches. The pipeline minimizes mismatches to inclusion targets while maximizing mismatches to non-targets. The algorithms were developed for production of Applied Biosystems' large collection of TaqMan® assays and have been extensively tested in-house to validate the assay design process.

The target selected by the assay design pipeline for *S. aureus* detection was determined to be unique in *S. aureus*. The assay was validated against an extensive inclusion and exclusion panel. The *S. aureus* assay is multiplexed with an internal positive control (IPC) which monitors PCR inhibition in each reaction.

The assay was tested in lyophilized and tableted format. Dried-down formats consisted of PCR reagent components, primers, and probes. *Staphylococcus* species were enriched overnight, lysed, and lysates added to the tablets in a 30 µL reaction. The reactions were run on DuPont Qualicon's BAX Q7 system which produced plus/minus calls and generated Ct (cycle threshold) values.




Figure 1. TaqMan® assay design pipeline

Specificity

| Inclusion | Isolates | Pos | Neg |
|-------------------------------------|----------|-----|-----|
| <i>Staphylococcus aureus</i> | 27 | 27 | 0 |
| Exclusion | | | |
| <i>Staphylococcus arletiae</i> | 1 | 0 | 1 |
| <i>Staphylococcus aureofaciens</i> | 1 | 0 | 1 |
| <i>Staphylococcus capitis</i> | 1 | 0 | 1 |
| <i>Staphylococcus caprae</i> | 2 | 0 | 2 |
| <i>Staphylococcus carnosus</i> | 4 | 0 | 4 |
| <i>Staphylococcus coelicolicus</i> | 1 | 0 | 1 |
| <i>Staphylococcus chromogenes</i> | 1 | 0 | 1 |
| <i>Staphylococcus cohnii</i> | 1 | 0 | 1 |
| <i>Staphylococcus delphini</i> | 1 | 0 | 1 |
| <i>Staphylococcus epidermidis</i> | 1 | 0 | 1 |
| <i>Staphylococcus felis</i> | 1 | 0 | 1 |
| <i>Staphylococcus gallinarum</i> | 1 | 0 | 1 |
| <i>Staphylococcus haemolyticus</i> | 1 | 0 | 1 |
| <i>Staphylococcus hominis</i> | 1 | 0 | 1 |
| <i>Staphylococcus intermedius</i> | 1 | 0 | 1 |
| <i>Staphylococcus saprophyticus</i> | 1 | 0 | 1 |
| <i>Staphylococcus schweitzeri</i> | 1 | 0 | 1 |
| <i>Staphylococcus sciuri</i> | 1 | 0 | 1 |
| <i>Staphylococcus simulans</i> | 1 | 0 | 1 |
| <i>Staphylococcus warneri</i> | 1 | 0 | 1 |
| <i>Staphylococcus xylosum</i> | 1 | 0 | 1 |

Table 1. The *S. aureus* assay was specific for all isolates of *S. aureus* in the inclusion panel tested. This included the detection of an isolate of a methicillin-resistant strain of *S. aureus*. None of the closely related *Staphylococcus* species in the exclusion panel were detected.






Figure 2. Amplification plots from *S. aureus* and non-*S. aureus* samples from the inclusion and exclusion panel.

Efficiency



Figures 3 and 4. Assay efficiency was determined by testing the assay over a dilution series of purified *S. aureus* genomic DNA. Each ten-fold DNA dilution was tested in quadruplicate. Efficiency was measured using the Ct slope method (plot of average Ct versus log DNA concentration). The expected slope for a ten-fold dilution series of DNA is -3.32 when efficiency is equal to 100% and would indicate a doubling of PCR product is occurring at every PCR cycle. The *S. aureus* assay demonstrated efficiency close to 100%.

Sensitivity

| <i>S. aureus</i> gDNA (pg) | Genomic copy number | Av Ct |
|----------------------------|---------------------|-------|
| 3000 | 1,000,000 | 22.58 |
| 300 | 100,000 | 25.83 |
| 30 | 10,000 | 29.40 |
| 3 | 1,000 | 32.83 |
| 0.3 | 100 | 35.70 |
| 0.03 | 10 | 37.43 |




Table 2 and Figure 5. Ten copies of *S. aureus* genomic DNA were reproducibly detected.

Comparison of Real-Time PCR to ISO reference method




Figure 6 and Table 3. Sensitivities of Real-Time PCR and the ISO reference method for detection of *S. aureus* were compared. Equivalent test samples were taken from twenty different 10 g portions of *S. aureus*-contaminated ground beef mixed with pre-washed Buttersfield buffer. The 1 mL test samples were separately (1) enriched in BHI/NaCl followed by lysis and Real-Time PCR and (2) plated on Baird Parker plates. Both Real-Time PCR and plating results were confirmed by latex agglutination and coagulase tests. The 1 mL test sample used for each method was estimated to contain on average 3 cfu *S. aureus*. Due to sampling distribution, it is expected that some samples would contain 0 cfu and would result in a negative result. The data in Table 3 show that the sensitivities of Real-Time PCR (18/20) and the ISO reference method (19/20) were comparable. Estimated concentration of the master sample were based on the test samples and are indicated in the final column of Table 3.

| Sample | Real-time PCR | | | | | ISO | | | | |
|--------|---------------|-----|---------|------------|-----------|-----|------------|-----------|-------|--|
| | BAX | BP | streaks | latex aggl | coagulase | BP | latex aggl | coagulase | cfu/g | |
| 1 | POS | POS | POS | POS | POS | 4+ | 4+ | 4+ | 40 | |
| 2 | POS | POS | POS | POS | POS | 5+ | 5+ | 5+ | 50 | |
| 3 | POS | POS | POS | POS | POS | 6+ | 6+ | 6+ | 60 | |
| 4 | POS | POS | POS | POS | POS | 4+ | 4+ | 4+ | 40 | |
| 5 | POS | POS | POS | POS | POS | 3+ | 3+ | 3+ | 30 | |
| 6 | POS | POS | POS | POS | POS | 4+ | 4+ | 4+ | 40 | |
| 7 | neg | neg | neg | neg | neg | 5+ | 5+ | 5+ | 50 | |
| 8 | POS | POS | POS | POS | POS | 3+ | 1+ | 1+ | 10 | |
| 9 | POS | POS | POS | POS | POS | 4+ | 3+ | 3+ | 30 | |
| 10 | POS | POS | POS | POS | POS | 4+ | 2+ | 2+ | 20 | |
| 11 | POS | POS | POS | POS | POS | 4+ | 2+ | 1+ | 10 | |
| 12 | POS | POS | POS | POS | POS | 4+ | 3+ | 3+ | 30 | |
| 13 | POS | POS | POS | POS | POS | 4+ | 4+ | 3+ | 30 | |
| 14 | POS | POS | POS | POS | POS | 3+ | 3+ | 1+ | 10 | |
| 15 | POS | POS | POS | POS | POS | 0+ | 0+ | 0+ | | |
| 16 | POS | POS | POS | POS | POS | 2+ | 0+ | 0+ | | |
| 17 | POS | POS | POS | POS | POS | 3+ | 3+ | 3+ | 30 | |
| 18 | POS | POS | POS | POS | POS | 4+ | 4+ | 4+ | 40 | |
| 19 | POS | POS | POS | POS | POS | 5+ | 1+ | 1+ | 10 | |
| 20 | neg | neg | neg | neg | neg | 3+ | 3+ | 3+ | 30 | |

Conclusions

- A highly specific TaqMan® assay for the detection of *S. aureus* was developed using Applied Biosystems' assay design pipeline.
- The sensitivity of the assay permits reliable detection of 10 genomic copies of *S. aureus*.
- The sensitivity of detection of *S. aureus* in contaminated ground beef is comparable to the sensitivity of the ISO reference method.
- The time to result—including enrichment, sample prep, and Real-Time PCR—is within 24 hrs.
- The available tableted form enables ease-of-use for rapid screening.

Trademark/Licensing

For Research Use Only. Not for use in diagnostic procedures.
 Practice of the patented 5' Nuclease Process requires a license from Applied Biosystems. The purchase of the TaqMan® assay includes an immunity from suit under patents specified in the product insert to use only the amount purchased for the purchaser's own internal research when used with the separate purchase of an Authorized 5' Nuclease Core Kit. No other patent rights are conveyed expressly, by implication, or by estoppel. For further information on purchasing licenses contact the Director of Licensing, Applied Biosystems, 850 Lincoln Center Drive, Foster City, California 94024, USA.
 Applied Biosystems, and AB (Design) are registered trademarks of Applied Biosystems Corporation or its subsidiaries in the US and/or other countries. BAX is a registered trademark of Qualicon, Inc. TaqMan is a registered trademark of Roche Molecular Systems, Inc. All other trademarks are the sole property of their respective owners.
 © 2007 Applied Biosystems. All rights reserved.