

GeneMapper® Software 5

Product	GeneMapper® Software
Version	5
OS Compatibility	Windows® 7 Professional, 32-bit edition, service pack 1
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	3130 series Data Collection 4 3730 series Data Collection 4 3500 Data Collection v.2
Stand Alone Computer Specs	Minimum: 2Ghz CPU or faster 2 Gb RAM 1 Gb Hard disk space for Client 7.5 Gb Hard Disk space for Full (500 Mb on boot drive, 7 Gb on drive where GeneMapper® Software is installed) Recommended: CPU: Core 2 Duo E8400/3.0Ghz 4 Gb RAM 1 Gb Hard disk space for Client 7.5 Gb Hard Disk space for Full (500 Mb on boot drive, 7 Gb on drive where GeneMapper® Software is installed)
File Type	.fsa
Sample File Compatibility	310 PC, Data Collection v.3.0, 3.1 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0, 3.1, 4 3730 series, Data Collection v. 3.0, 3.1, 4 3500 series, Data Collection v.1.0, v.1.1, 2

Notes

This software has been extensively tested on Windows® 7 SP1 (32-bit) version. While there are some indications that software will run normally on the 64-bit version of Windows® 7 SP1 operating system, the 64-bit configuration is not officially supported. It is strongly recommended that this software is run on Windows® 7 SP1 (32-bit). If installing on the instrument computer, the update will require a new computer and upgrading the Data Collection Software as well.