

# GeneMapper® Software v.4.1

Product	GeneMapper® Software
Version	4.1
OS Compatibility	Windows® XP Professional, service packs 2 and 3 Windows® Vista® (32-bit Business Edition), service pack 1
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	3130 series Data Collection v.3.0 3730 series Data Collection v.3.0 310 Data Collection v.3.1
Stand Alone Computer Specs	Minimum: 2.2Ghz CPU or faster 1 Gb RAM 160 Gb Hard disk space * Recommended: CUP: Core 2 Duo E8400/3.0 Ghz 4 Gb RAM 320 Gb Hard Disk Space*
File Type	.fsa
Sample File Compatibility	310 PC, Data Collection v.3.0, 3.1 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0, 3.1, 4 3730 series, Data Collection v. 3.0, 3.1, 4 3500 series, Data Collection v.1.0, v.1.1, 2
Notes	Upgrading the Data Collection software on the 3130 and 3730 series instruments to v.3.1 for co-installation is required. Software is not Windows® 7 compatible. * Hard drive requirements will vary depending on how much space is allocated to the Oracle Database. Having a large hard drive will allow for database expansion should the database fill up.