

# GeneMapper® Software v.4.0

Product	GeneMapper® Software
Version	4.0
OS Compatibility	Windows® XP Professional, service pack 1 or 2 Windows® 2000, service pack 4
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	3130 series Data Collection v.3.0 3730 series Data Collection v.3.0 310 Data Collection v.3.1
Stand Alone Computer Specs	Minimum: 2.4Ghz CPU or faster 1 Gb RAM 120 Gb Hard disk space *
File Type	.fsa
Sample File Compatibility	310 PC, Data Collection v.3.0, 3.1 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0, 3.1, 4 3730 series, Data Collection v. 3.0, 3.1, 4
Notes	Software is not Windows® 7 compatible. * Hard drive requirements will vary depending on how much space is allocated to the Oracle Database. Having a large hard drive will allow for database expansion should the database fill up.