

SeqScape® Software v.2.7

Product	SeqScape® Software
Version	2.7
OS Compatibility	Windows® XP Professional, service pack 2 or 3 Windows® Vista® (32-bit Business Edition), service pack 1
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	310 PC, Data Collection v.3.1 3130 series Data Collection v.3.1 3730 series Data Collection v.3.1 3500 series Data Collection v.1.0
Stand Alone Computer Specs	Minimum: 2.8 Ghz Pentium® IV CPU or faster 1 Gb RAM 40 Gb Hard disk space for application Recommended: Core 2 Duo E8400/3.0Ghz CPU 2 Gb RAM 320 Gb Hard disk space
File Type	.ab1
Sample File Compatibility	310 PC, Data Collection v.3.0, 3.1 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0, 3.1, 4 3730 series, Data Collection v. 3.0, 3.1, 4 3500 series, Data Collection v.1.0, v.1.1, 2
Notes	*Upgrade will require Data Collection Software v.3.1 on the 3130 and 3730 series instruments for co- installation. System upgrade may be required. *Software is not Windows® 7 compatible.

Notes

*Demo version is no longer available through Life Technologies™, only through third party sites. Life Technologies™ cannot guarantee the integrity of the software or that it is virus-free on these sites, use at own risk.