

SeqScape® Software v.3

Product	SeqScape® Software
Version	3
OS Compatibility	Windows® 7 Professional, 32-bit edition, service pack 1
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	3130 series Data Collection v.4 3730 series Data Collection v.4 3500 series Data Collection v.2
Stand Alone Computer Specs	Minimum: 2.0 Ghz or faster Intel processor 2 Gb RAM 1 Gb Hard disk space for application Recommended: Core 2 Duo E8400/3.0Ghz CPU 4 Gb RAM 1 Gb Hard disk space
File Type	.ab1
Sample File Compatibility	310 PC, Data Collection v.3.0, 3.1 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0, 3.1, 4 3730 series, Data Collection v. 3.0, 3.1, 4 3500 series, Data Collection v.1.0, v.1.1, 2
Notes	This software has been extensively tested on Windows® 7 SP1 (32-bit) version. While there are some indications that the software will run normally on the 64-bit version of Windows® 7 SP1 operating system, the 64-bit configuration is not officially

Notes

supported. It is strongly recommended that this software is run on Windows® 7 SP1 (32-bit). If installing on the instrument computer, the update will require a new computer and upgrading the Data Collection Software as well.