

Sequencing Analysis Software 6

Product	Sequencing Analysis
Version	6
OS Compatibility	Windows® 7 Professional, 32-bit edition, service pack 1
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	3130 series Data Collection v.4 3730 series Data Collection v.4 3500 Data Collection v.2
Stand Alone Computer Specs	Minimum: 2Ghz CPU or faster 2 Gb RAM 1 Gb Hard disk space for application Recommended: Core 2 Duo E8400/3.0Ghz CPU 4 Gb RAM 1 Gb Hard disk space
File Type	.ab1
Sample File Compatibility	310 PC, Data Collection v.3.0, 3.1 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0, 3.1, 4 3730 series, Data Collection v. 3.0, 3.1, 4 3500 series, Data Collection v.1.0, v.1.1, 2
Notes	This software has been extensively tested on Windows® 7 SP1 (32-bit) version. While there are some indications that the software will run normally on the 64-bit version of Windows® 7 SP1 operating system, the 64-bit configuration is not officially supported.

Notes

It is strongly recommended that this software is run on Windows® 7 SP1 (32-bit).

If installing on the instrument computer, the update will require a new computer and upgrading the Data Collection Software as well.