

Sequencing Analysis Software v.5.4

Product	Sequencing Analysis
Version	5.4
OS Compatibility	Windows® XP Professional, service pack 2 or 3 Windows® Vista® (32-bit Business Edition), service pack 1
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	3130 series Data Collection v.3.1 3730 series Data Collection v.3.1 310 Data Collection v.3.1 3500 series Data Collection 1.0
Stand Alone Computer Specs	Minimum: 2.8 Ghz Pentium® IV CPU or faster 1 Gb RAM 40 Gb Hard disk space for application Recommended: Core 2 Duo E8400/3.0Ghz CPU 2 Gb RAM 320 Gb Hard disk space
File Type	.ab1
Sample File Compatibility	310 PC, Data Collection v.3.0, 3.1 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0, 3.1, 4 3730 series, Data Collection v. 3.0, 3.1, 4 3500 series, Data Collection v.1.0, v1.1, 2
Notes	If installing on the instrument computer, the upgrade will require upgrading to Data Collection Software v.3.1 on the 3130 and 3730 series instruments. Software is not Windows® 7 compatible.