

# Sequencing Analysis Software v.5.3.1

Product	Sequencing Analysis
Version	5.3.1
OS Compatibility	Windows® 2000, service pack 3 or 4 Windows® XP Professional, service pack 2 or 3 Windows® Vista® (32-bit Business Edition)
Part Numbers	Please contact your local Sales Representative for more information.
Instrument Installation Compatibility	3130 series Data Collection v.3.0 3730 series Data Collection v.3.0 310 Data Collection v.3.1
Stand Alone Computer Specs	Minimum: 733 Mhz Pentium® IV CPU or faster 512 Mb RAM 5 Gb Hard disk space for application
File Type	.ab1
Sample File Compatibility	310 PC, Data Collection v.3.0 3100 series, Data Collection v.1.0, 1.1, 2.0 3130 series, Data Collection v. 3.0 3730 series, Data Collection v. 3.0
Notes	Software is not Windows® 7 compatible.