

XL Core Version 1.0.115 open source software components

This document describes the source code used on the XL Core with Version 1.0.115 of the embedded software, which falls under the GNU public license (the GPL) or the GNU Lesser General Public License (the LGPL).

For more information about the GPL, see: <http://www.gnu.org/licenses/gpl-2.0.html>

For more information about the LGPL, see: <http://www.gnu.org/licenses/lgpl-2.1.html>

Source Code download

Package	License	Homepage
Linux Kernel 2.6.37 with AMG patches	GPLv2	http://www.kernel.org
Linux Kernel Config 2.6.37	GPLv2	http://www.kernel.org

Toolchain

The entire XL Core Software for version 1.0.115 was built using the following Code Sourcery tools and libc:

Sourcery G++ Lite 2009q1-203 for ARM GNU/Linux available [here](#), itself based on the following components:

- GNU binutils-2.19.51
- gcc-4.3.3
- eglibc-2.8
- gdb-6.8.50