


Nunc Flasks, Dishes & Multidishes

Buy any 2 cases of qualified Nunc cell culture-treated products and receive 1 case of your choice at no additional cost!

Achieve consistent and reproducible cell growth

Thermo Scientific Nunc Flasks, Dishes and Multiwell Plates are treated with the proprietary Thermo Scientific™ Nunclon™ Delta surface that is designed to provide maximum adhesion and growth for a broad range of cell types.

- Thermo Scientific™ Nunc™ EasYFlasks™ feature an angled neck and wide opening, allowing easy access to cells by pipettes or cell scrapers
- Nunc Dishes and Multidishes are treated for adherent cultures and deliver optimal clarity, quality and sterility. Multiwell plates feature raised well rims that help to mitigate contamination risks



Thermo Scientific Nunclon Delta treated flasks, dishes or multidishes

Purchase 2 cases of the following Nunclon flasks, dishes & multidishes

Cat. No.	Description	No./ Pack/ Case
Thermo Scientific Nunc EasyFlasks		
156367	EasyFlask, Nunclon Delta-treated, Filter cap, 25 cm ² , 7 mL working volume	10/200
156499	EasyFlask, Nunclon Delta-treated, Filter cap, 75 cm ² , 25 mL working volume	5/100
159910	EasyFlask, Nunclon Delta-treated, Filter cap, 175 cm ² , 55 mL working volume	5/30
159934	EasyFlask, Nunclon Delta-treated, Filter cap, 225 cm ² , 70 mL working volume	5/30
156340	EasyFlask, Nunclon Delta-treated, Vent/close cap, 25 cm ² , 7 mL working volume	10/200
156472	EasyFlask, Nunclon Delta-treated, Vent/close cap, 75 cm ² , 25 mL working volume	5/100
159920	EasyFlask, Nunclon Delta-treated, Vent/close cap, 175 cm ² , 55 mL working volume	5/30
159933	EasyFlask, Nunclon Delta-treated, Vent/close cap, 225 cm ² , 70 mL working volume	5/30
Thermo Scientific Nunc TripleFlasks		
132867	TripleFlask, Nunclon Delta-treated, Vent/close cap, 500 cm ² , 200 mL working volume	2/32
132913	TripleFlask, Nunclon Delta-treated, Filter cap, 500 cm ² , 200 mL working volume	2/32
132920	TripleFlask with Bar Code, Nunclon Delta-treated, Filter cap, 500 cm ² , 200 mL working volume	2/32

Additional skus on the next page >


Redeem now at
thermofisher.com/cellculturepromos

Terms and Conditions:

This promotion is open to customers in the U.S. and Canada only. Valid for orders placed between January 1, 2017 and June 30, 2017. Redemptions for the promotional case and proof of purchase must be submitted by June 30, 2017. Promotion offers for storage plates and tubes are limited to 5 redemptions per laboratory. Subject to dealer participation. Cannot be combined with other discounts or promotion. Offer void where prohibited, licensed or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply. U.S. customers acknowledge that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

For Research Use Only. Not for use in diagnostic procedures. © 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. Cannot be combined with other discounts or promotion. **LSPCCH1Promo0117**

Nunc Flasks, Dishes & Multidishes



Nunclon Delta

Unique to our Thermo Scientific Nunc cell culture plastic is our **Nunclon Delta surface** to promote maximum adhesion for a broad range of cell types.

The Nunclon Delta surface is a fully synthetic energy-treated surface that makes the otherwise very hydrophobic polystyrene surface more hydrophilic, thus facilitating cell attachment and growth. The Nunclon Delta certificate ensures high quality standards and lot-to-lot consistency.

(cont'd from page 4)

Purchase 2 cases of the following flasks, dishes & multidishes		
Cat. No.	Description	No./ Pack/ Case
Thermo Scientific Nunc		
153066	35 mm Dish, Nunclon Delta-treated, 8.8 cm ² , with lid; with airvent	10/500
150288	60 mm Dish, Nunclon Delta-treated, 21.5 cm ² , with lid; with airvent	10/400
172931	100x20mm Dish, Nunclon Delta-treated, 56.7 cm ² , with lid; with airvent	10/240
Thermo Scientific Nunc Dishes		
168381	150x20mm Dish, Nunclon Delta-treated, 145 cm ² , with lid; with airvent	10/80
150350	100X15 mm Dish, Nunclon Delta-treated, 56.7 cm ² , with lid; with airvent	150
157150	150x15 mm Dish, Nunclon Delta-treated, 145 cm ² , with lid; airvent	10/120
140675	6-well Multidish, Nunclon Delta-treated, 9.6 cm ² , round well; with lid	1/75
140685	6-well Multidish, Nunclon Delta-treated, 9.6 cm ² , round well; with lid	5/85
150628	12-well Multidish, Nunclon Delta-treated, 3.5cm ² , round well; with lid; sterile	1/75
142475	24-well Multidish, Nunclon Delta-treated, 1.8 cm ² , round well; with lid; sterile	1/75
142485	24-well Multidish, Nunclon Delta-treated, 1.8 cm ² , round well; with lid; sterile	5/85
150687	48-well Multidish, Nunclon Delta-treated, 1.1 cm ² , round well; with lid; sterile	1/75
152640	48-well Multidish, Nunclon Delta-treated, 1.1 cm ² , round well; with lid; sterile	5/85
Thermo Scientific Nunc Dishes		
172931	100 mm Dish, Nunclon Delta-treated, 56.7 cm ² , with lid; with airvent	10/240
150350	100 mm Dish, Nunclon Delta-treated, 56.7 cm ² , with lid; with airvent	150/150
168381	150 mm Dish, Nunclon Delta-treated, 145 cm ² , with lid; with airvent	10/80
157150	150 mm Dish, Nunclon Delta-treated, 150 cm ² , with lid; with airvent	10/120
166508	245 mm Dish, Nunclon Delta-treated, 500 cm ² , with lid; without airvent	4/16
Thermo Scientific Nunc Multidishes		
176740	4-well Multidish, Nunclon Delta-treated, 1.9 cm ² , round well; with lid	4/120
167063	4-well Multidish, Nunclon Delta-treated, 21.8 cm ² , round well; with lid	10/100
140675	6-well Multidish, Nunclon Delta-treated, 9.6 cm ² , round well; with lid	1/75
140685	6-well Multidish, Nunclon Delta-treated, 9.6 cm ² , round well; with lid	5/85
167064	8-well Multidish, Nunclon Delta-treated, 10.5 cm ² , rectangular well; with lid; sterile	10/100
150628	12-well Multidish, Nunclon Delta-treated, 3.5 cm ² , round well; with lid; sterile	1/75
142475	24-well Multidish, Nunclon Delta-treated, 1.8 cm ² , round well; with lid; sterile	1/75
142485	24-well Multidish, Nunclon Delta-treated, 1.8 cm ² , round well; with lid; sterile	5/85
150687	48-well Multidish, Nunclon Delta-treated, 1.1 cm ² , round well; with lid; sterile	1/75
152640	48-well Multidish, Nunclon Delta-treated, 1.1 cm ² , round well; with lid; sterile	5/85

Terms and Conditions:

This promotion is open to customers in the U.S. and Canada only. Valid for orders placed between January 1, 2017 and June 30, 2017. Redemptions for the promotional case and proof of purchase must be submitted by June 30, 2017. Promotion offers for storage plates and tubes are limited to 5 redemptions per laboratory. Subject to dealer participation. Cannot be combined with other discounts or promotion. Offer void where prohibited, licensed or restricted by federal, state, provincial, or local laws or regulation or agency/institutional policy. Other restrictions may apply. U.S. customers acknowledge that this offer may include a discount or other price reduction that must be properly and accurately accounted for and reported by customer in accordance with all federal and state laws, including without limitation the federal anti-kickback law (42 C.F.R. § 1320a-7b(b)(3)(A)) and regulations thereunder (42 C.F.R. §1001.952(h)).

For Research Use Only. Not for use in diagnostic procedures. © 2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. Cannot be combined with other discounts or promotion. **LSPCCH1Promo0117**