



Rat IgG1 kappa Isotype Control (eBRG1), Super Bright™ 702, eBioscience™

Product Details		
Size	100 μg	
Host/Isotype	Rat / IgG1, kappa	
Class	Monoclonal	
Туре	Isotype Control	
Clone	eBRG1	
Conjugate	Super Bright™ 702	
Excitation/Emission Max	413/702 nm	
Form	Liquid	
Concentration	0.2 mg/mL	
Purification	Affinity chromatography	
Storage buffer	PBS, pH 7.2, with BSA	
Contains	0.09% sodium azide	
Storage conditions	4° C, store in dark, DO NOT FREEZE!	
RRID	AB_2717169	

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	Assay-Dependent	-
Control (Ctrl)	Assay-Dependent	-

Product Specific Information

Description: The monoclonal rat IgG1, kappa is useful as an isotype control immunoglobulin.


Applications Reported: This rat IgG1 isotype control has been reported for use in flow cytometric analysis and ELISA.


Applications Tested: This rat IgG1 isotype control has been used by flow cytometric analysis of mouse splenocyte suspension, This can be used at the same concentration as the experimental antibody.

Excitation: 405 nm; Emission: 702 nm; Laser: Violet Laser.

When using two or more Super Bright dye-conjugated antibodies in a staining panel, it is recommended to use Super Bright Complete Staining Buffer (Product # SB-4401) to minimize any non-specific polymer interactions. Please refer to the datasheet for Super Bright Staining Buffer for more information.

Product Images For Rat IgG1 kappa Isotype Control (eBRG1), Super Bright™ 702, eBioscience™





Rat IgG1 kappa Isotype Control (67-4301-82)

Staining of C57Bl/6 splenocytes with Anti-Mouse CD4 FITC (Product # 11-0041-82) and staining buffer (autofluorescence, left) or 0.06 µg of Rat IgG1 kappa Isotype Control Super Bright 702 (right). Total viable cells were used for analysis.

□1 Reference

An engineered IL-2 partial agonist promotes CD8+ T cell stemness. Nature (2021)

Propofol attenuates inflammatory damage on neurons following cerebral infarction by inhibiting excessive activation of microglia. Int J Mol Med (2019)

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production obcumentation, specifications and/or accompanying package inserts ("Documentation,"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is imited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample increase in series illustrated the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT.

BUYER'S EXCLUSIVE REMEDY FOR NON-CORNOR/MINO PRODUCTS OURING THE WARRANTY PERFIOI DIS LIMITACE COR REFUND FOR THE HE NON-CONFORMING PRODUCTS; AT SELLER'S SOLE OPTION, THERE IS NO OBLICATION TO REPAIR, REPLACE OR REFUND FOR THE HE NON-CONFORMING PRODUCTS; AT SELLER'S SOLE OPTION, THERE IS NO OBLICATION TO REPAIR, REPLACE OR REFUND FOR THE PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR