CD4 Monoclonal Antibody (GK1.5), PE-Cyanine5, eBioscience™

Product Details

i loudet Details	
Size	100 µg
Species Reactivity	Mouse
Published Species	Fungi, Rat, Fish, Mouse, Human
Host/Isotype	Rat / IgG2b, kappa
Recommended Isotype Control	Rat IgG2b kappa Isotype Control (eB149/10H5), PE-Cyanine5, eBioscience™
Class	Monoclonal
Туре	Antibody
Clone	GK1.5
Conjugate	PE-Cyanine5
Excitation/Emission Max	568/666 nm
Form	Liquid
Concentration	0.2 mg/mL
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.2
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_468695

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	0.06 μg/test	197 Publications
Functional Assay (FN)	-	2 Publications

Product Specific Information

Description: The GK1.5 monoclonal antibody reacts with the mouse CD4 molecule, a 55 kDa cell surface receptor expressed by a majority of thymocytes, subpopulation of mature T cells and dendritic cells. CD4 binds to MHC class II on the surface of antigen presenting cells and plays an important role both in T cell development and in optimal functioning of mature T cells. In T cells, CD4 associates with protein tyrosine kinase p56lck through its cytoplasmic tail. Binding of GK1.5 is blocked by RM4-5.

Applications Reported: The GK1.5 antibody has been reported for use in flow cytometric analysis.

Applications Tested: The GK1.5 antibody has been tested by flow cytometric analysis of mouse thymocytes and splenocytes. This can be used at less than or equal to 0.06 μ g per test. A test is defined as the amount (μ g) of antibody that will stain a cell sample in a final volume of 100 μ L. Cell number should be determined empirically but can range from 10^5 to 10^8 cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.

Light sensitivity: This tandem dye is sensitive photo-induced oxidation. Please protect this vial and stained samples from light.


Fixation: Samples can be stored in IC Fixation Buffer (cat. 00-8222) (100 μ L cell sample + 100 μ L IC Fixation Buffer) or 1-step Fix/Lyse Solution (cat. 00-5333) for up to 3 days in the dark at 4°C with minimal impact on brightness and FRET efficiency /compensation. Some generalizations regarding fluorophore performance after fixation can be made, but clone specific performance should be determined empirically.

1

Excitation: 488-561 nm; Emission: 667 nm; Laser: Blue Laser, Green Laser, Yellow-Green Laser.


Filtration: 0.2 µm post-manufacturing filtered.

Product Images For CD4 Monoclonal Antibody (GK1.5), PE-Cyanine5, eBioscience™



CD4 Antibody (15-0041-82) in Flow

Staining of C57Bl/6 splenocytes with staining buffer (autofluorescence) (open histogram) or 0.03 μ g of Anti-Mouse CD4 PE-Cyanine5 (filled histogram).Total cells were used for analysis.



CD4 Antibody (15-0041-82)

Staining of mouse splenocytes. As expected based on known relative expression patterns, CD4 clone GK1.5 stains about 25% of total lymphocytes and does not stain CD45R (B220)+ B cells. Details: Balb/c splenocytes were surface stained with CD45R (B220, clone RA3-6B2) and Rat IgG2b Isotype Control or CD4 (clone GK1.5). Cells in the lymphocyte gate were used for analysis. {RE}

View more figures on thermofisher.com

□ 199 References

Flow Cytometry (197)

Theranostics		Year 2023
HSF1 promotes CD69 ⁺ Treg differentiation to inhibit colitis progression.		
"Published figure using CD4 monoclonal antibody (Product # 15-0041-82) in Flow Cytometry"		
Authors: Yu L,Zhou B,Zhu Y,Li L,Zhong Y,Zhu L,Wang H,Chen H,Xu J,Guo T,Feng L,Wang X,Cai Z,Wang J,Jin H		
Journal of immunology (Baltimore, Md. : 1950)		Year
Prolidase Deficiency Causes Spontaneous T Cell Activation and Lupus-		2023
like Autoimmunity.		
"Published figure using CD4 monoclonal antibody (Product # 15-0041-82) in Flow Cytometry"		
Authors: Hodgson R,Crockford TL,Bhandari A,Kepple JD,Back J,Cawthorne E,Abeler-Dörner L,Laing AG,Clare S, Speak A,Adams DJ,Dougan G,Hayday AC,Deobagkar-Lele M,Cornall RJ,Bull KR		
	View more Flow references	on thermofisher.co
unctional Assay (2)		
Investigative ophthalmology & visual science		Year
Suppression of alkali burn-induced corneal neovascularization by		2008
dendritic cell vaccination targeting VEGF receptor 2.		

Authors: Mochimaru H,Usui T,Yaguchi T,Nagahama Y,Hasegawa G,Usui Y,Shimmura S,Tsubota K,Amano S, Kawakami Y,Ishida S

Circulation Myocardial infarct-sparing effect of adenosine A2A receptor activation is due to its action on CD4+ T lymphocytes. Authors: Yang Z,Day YJ,Toufektsian MC,Xu Y,Ramos SI,Marshall MA,French BA,Linden J	Year 2006
--	---------------------

More applications with references on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individual. Unless otherwise stated in the Documentation, this warranty for on or even from date of shipment when the Product is subject to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample inneally likes to the general type and quality of goods and does not represent that any Product will conform to such model or sample inneaded usage. This warranty does not extend to anyone other than the Buyer. Any model or sample inneally likes to the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTY EXICULAR PURPOSE, OR NON INFRINGEMENT. BUYERYS EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS ON UNING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCTS) ASTER OR EVENT OF FOR EVENT OF FOR REFUND FOR THE NON-CONFORMING PRODUCTS IN THERE IS NOL OPTION. THERE IS NO OBLIGATION TO FEPAIR, REPLACE OR REFUND FOR THE RONL-CONFORMING PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HADLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.