

CD117 (c-Kit) Monoclonal Antibody (2B8), APC-eFluor™ 780, eBioscience™

Product Details	
Size	100 µg
Species Reactivity	Mouse, Pig
Published Species	Human, Mouse
Host/Isotype	Rat / IgG2b, kappa
Recommended Isotype Control	Rat IgG2b kappa Isotype Control (eB149/10H5), APC-eFluor™ 780, eBioscience™
Class	Monoclonal
Type	Antibody
Clone	2B8
Conjugate	APC-eFluor™ 780
Excitation/Emission Max	756/785 nm
Form	Liquid
Concentration	0.2 mg/mL
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.2
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_1272177

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	0.125 µg/test	86 Publications
Functional Assay (FN)	-	1 Publication
In vitro Assay (IV)	-	1 Publication

Product Specific Information

Description: The 2B8 monoclonal antibody reacts with mouse CD117, also known as c-Kit receptor, Steel factor receptor, and stem cell factor receptor. A member of the tyrosine kinase receptor family, this 145-kDa molecule is expressed by a majority of hematopoietic progenitor cells characterized in the mouse bone marrow as a small subset of cells positive for Sca-1 and Thy1 (Thy1 low) and negative for lineage markers. The interaction of the mouse c-Kit receptor and steel factor promotes the proliferation and differentiation of hematopoietic progenitor cells. CD117 is also expressed by mast cells and plays a role in signaling and activation of these cells.

Applications Reported: This 2B8 antibody has been reported for use in flow cytometric analysis.

Applications Tested: This 2B8 antibody has been tested by flow cytometric analysis of mouse bone marrow cells. This can be used at less than or equal to 0.125 µg per test. A test is defined as the amount (µg) of antibody that will stain a cell sample in a final volume of 100 µL. Cell number should be determined empirically but can range from 10⁵ to 10⁸ cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.

APC-eFluor 780 emits at 780 nm and is excited with the Red laser (633 nm). Please make sure that your instrument is capable of detecting this fluorochrome.


Light sensitivity: This tandem is sensitive to photo-induced oxidation. Please protect this vial and stained samples from light.

Fixation: Samples can be stored in IC Fixation Buffer (cat. 00-8222) (100 µL cell sample + 100 µL IC Fixation Buffer) or 1-step Fix/Lyse Solution (cat. 00-5333) for up to 3 days in the dark at 4°C with minimal impact on brightness and FRET efficiency /compensation. Some generalizations regarding fluorophore performance after fixation can be made, but clone specific performance should be determined empirically.

Excitation: 633-647 nm; Emission: 780 nm; Laser: Red Laser.

Filtration: 0.2 µm post-manufacturing filtered.

Product Images For CD117 (c-Kit) Monoclonal Antibody (2B8), APC-eFluor™ 780, eBioscience™



CD117 (c-Kit) Antibody (47-1171-82) in Flow
Staining of C57BL/6 bone marrow cells with Anti-Mouse Ly-6A/E (Sca-1) PE (Product # 12-5981-82) and staining buffer (autofluorescence) (left) or 0.06 µg of Anti-Mouse CD117 (c-Kit) APC-eFluor® 780 (right). Cells in the Lineage (low) gate were used for analysis.

Flow Cytometry (86)

<p>Frontiers in immunology</p> <p>Insulin-regulated aminopeptidase contributes to setting the intensity of FcR-mediated inflammation.</p> <p>"47-1171-82 was used in Flow cytometry/Cell sorting to illustrate that IRAP-positive endosomal compartments, in promoting SHP-1 inhibition during FcR signaling, control the extent of phosphorylation events at the plasma membrane and contribute to setting the intensity of immune-complex triggered inflammatory diseases."</p> <p>Authors: Bratti M,Vibhushan S,Longé C,Koumantou D,Ménasché G,Benhamou M,Varin-Blank N,Blank U,Saveanu L,Ben Mkaddem S</p>	<p>Year 2022</p> <p>Species Mouse</p>
<p>Cell reports</p> <p>PDGFR⁺ cells play a dual role as hematopoietic precursors and niche cells during mouse ontogeny.</p> <p>"47-1171-82 was used in Flow cytometry/Cell sorting to demonstrate that PDGFR⁺ cells play a dual role in murine hematopoiesis."</p> <p>Authors: Sá da Bandeira D,Kilpatrick AM,Marques M,Gomez-Salazar M,Ventura T,Gonzalez ZN,Stefancova D,Rossi F,Vermeren M,Vink CS,Beltran M,Henderson NC,Jung B,van der Linden R,van de Werken HJG,van Ijcken WFJ,Betsholtz C,Forbes SJ,Cuervo H,Crisan M</p>	<p>Year 2022</p> <p>Species Mouse</p> <p>Dilution 1:800</p>

[View more Flow references on thermofisher.com](#)

Functional Assay (1)

<p>eLife</p> <p>Mice deficient of <i>Myc</i> super-enhancer region reveal differential control mechanism between normal and pathological growth.</p> <p>"47-1171 was used in Functional assays to show that a deletion in the gene desert upstream of the MYC oncogene on chromosome 8q24 can confer tumourigenic resistance."</p> <p>Authors: Dave K,Sur I,Yan J,Zhang J,Kaasinen E,Zhong F,Blaas L,Li X,Kharazi S,Gustafsson C,De Paepe A,Månsson R,Taipale J</p>	<p>Year 2017</p> <p>Species Mouse</p>
---	---

In vitro Assay (1)

<p>Blood</p> <p>HIF prolyl hydroxylase 2 (PHD2) is a critical regulator of hematopoietic stem cell maintenance during steady-state and stress.</p> <p>"47-1171 was used in Flow cytometry/Cell sorting to study the role of PHD2 in HSC compartment maintenance."</p> <p>Authors: Singh RP,Franke K,Kalucka J,Mamlouk S,Muschter A,Gembarska A,Grinenko T,Willam C,Naumann R,Anastassiadis K,Stewart AF,Bornstein S,Chavakis T,Breier G,Waskow C,Wielockx B</p>	<p>Year 2013</p> <p>Species Mouse</p>
---	---

More applications with references on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON-INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.