

# CD197 (CCR7) Monoclonal Antibody (3D12), Super Bright™ 645, eBioscience™

Product Details	
Size	100 Tests
Species Reactivity	Human
Host/Isotype	Rat / IgG2a, kappa
Recommended Isotype Control	Rat IgG2a kappa Isotype Control (eBR2a), Super Bright™ 645, eBioscience™
Class	Monoclonal
Туре	Antibody
Clone	3D12
Conjugate	Super Bright™ 645
Excitation/Emission Max	414/645 nm
Form	Liquid
Concentration	5 μL/Test
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.2, with BSA
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_2744870

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	5 μL (1 μg)/test	5 Publications

## **Product Specific Information**

Description: The 3D12 monoclonal antibody reacts with human CCR7, also known as EBI-1 and CD197. CCR7 is a member of the G-protein-coupled chemokine receptor family with seven membrane-spanning domains and functions as a receptor for 6Ckine/SLC (secondary lymphoid-tissue chemokine), CCL19 and CCL21. CCR7 has been shown to be internalized via clathrincoated pits and the majority recycled back to the plasma membrane. CCR7 is expressed on T cells and can be used to distinguish populations of naive from central and effector memory T cells. CCR7 has been shown to play a role in migration of memory T cells to inflamed tissue. Expression of CCR7 is also found on DC's. During DC maturation CCR7 expression increases and is thought to be involved in a variety of functions: chemotaxis to the lymph node, cellular architecture, rate of endocytosis, survival and maturation. Expression of CCR7 on the cell surface can be down regulated upon ligand binding.

It is recommended that the staining incubation time be increased to at least 45 minutes at 2-8°C for optimal staining.

Applications Reported: This 3D12 antibody has been reported for use in flow cytometric analysis.

Applications Tested: This 3D12 antibody has been pre-diluted and tested by flow cytometric analysis of normal human peripheral blood cells. This may be used at 5 µL (1 µg) per test. A test is defined as the amount (µg) of antibody that will stain a cell sample in a final volume of 100 µL. Cell number should be determined empirically but can range from 10^5 to 10^8 cells

Super Bright 645 is a tandem dye that can be excited with the violet laser line (405 nm) and emits at 645 nm. We recommend using a 660/20 bandpass filter. Please make sure that your instrument is capable of detecting this fluorochrome.

When using two or more Super Bright dye-conjugated antibodies in a staining panel, it is recommended to use Super Bright Complete Staining Buffer (Product # SB-4401) to minimize any non-specific polymer interactions. Please refer to the datasheet for Super Bright Staining Buffer for more information.


Light sensitivity: This tandem dye is sensitive to photo-induced oxidation. Please protect this vial and stained samples from light.


Fixation: Samples can be stored in IC Fixation Buffer (Product # 00-8222) (100 µL of cell sample + 100 µL of IC Fixation Buffer) or 1-step Fix/Lyse Solution (Product # 00-5333) for up to 3 days in the dark at 4°C with minimal impact on brightness and FRET efficiency/compensation. Some generalizations regarding fluorophore performance after fixation can be made, but clone specific performance should be determined empirically.

Excitation: 405 nm; Emission: 645 nm; Laser: Violet Laser

Super Bright Polymer Dyes are sold under license from Becton, Dickinson and Company.

# Product Images For CD197 (CCR7) Monoclonal Antibody (3D12), Super Bright™ 645, eBioscience™





## CD197 (CCR7) Antibody (64-1979-42) in Flow

Normal human peripheral blood cells were stained with CD45RA Monoclonal Antibody, FITC (Product # 11-0458-42) and Rat IgG2a kappa Isotype Control, Super Bright 645 (Product # 64-4321-82) (left) or CD197 (CCR7) Monoclonal Antibody, Super Bright 645 (right). Cells in the lymphocyte gate were used for analysis.

## View more figures on thermofisher.com

### □ 5 References

### Flow Cytometry (5)

Journal of inflammation research

Alteration of the Immune Microenvironment in HBsAg and HBeAg Dual-

"Published figure using CD197 (CCR7) monoclonal antibody (Product # 64-1979-42) in Flow Cytometry" Authors: Gao F, Wang H, Li X, Guo F, Yuan Y, Wang X, Zhang Y, Bai G

Positive Pregnant Women Presenting a High HBV Viral Load.

Frontiers in immunology

Neuroimmune Consequences of elF4E Phosphorylation on Chemotherapy-Induced Peripheral Neuropathy.

"Published figure using CD197 (CCR7) monoclonal antibody (Product # 64-1979-42) in Flow Cytometry"

Authors: Agalave NM, Mody PH, Szabo-Pardi TA, Jeong HS, Burton MD

Year 2022

Year 2021

View more Flow references on thermofisher.com

### More applications with references on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package interest ("Documentation"). No claim of suitability for use in applications to sample, to small or sample, the product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. No OTHER WARRANTIES, EXPRESS OR IMPLED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES FOR ANY PARTICULAR PURPOSE, OR NON INFIRIORIEMENT.

BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY EPIGO IS LIMITED TO REPAIR, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS; SOLE DETION, THERE IS NO BOLIGATION TO REPAIR, REPLACE OR REFUND FOR THE PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR REGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, ur diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals