# TNF alpha Monoclonal Antibody (68B 3C5), Biotin

<b>Product Details</b>	
Size	100 μg
Species Reactivity	Human
Published Species	Dog, Human
Host/Isotype	Mouse / IgG1
Class	Monoclonal
Туре	Antibody
Clone	68B 3C5
Conjugate	Biotin
Immunogen	Purified recombinant human TNF-alpha.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Protein A
Storage buffer	PBS, pH 7.4, with 5mg/mL BSA
Contains	0.5mg/mL sodium azide
Storage conditions	4° C
RRID	AB_2536272

Applications	Tested Dilution	Publications
ELISA (ELISA)	Assay-dependent	12 Publications

## **Product Specific Information**

This product is recommended as a detection antibody in Sandwich ELISA applications.

This clone can also be referred to as clone number 68B 3C5 4B3 A6.

#### □ 12 References

## **ELISA (12)**

Communications biology

Analysis of the glyco-code in pancreatic ductal adenocarcinoma identifies glycan-mediated immune regulatory circuits.

"AHC3419 was used in Enzyme-linked immunosorbent assay to show mechanisms by which the glyco-code contributes to the tolerogenic microenvironment in PDAC."

Authors: Rodriguez E,Boelaars K,Brown K,Maduni K,van Ee T,Dijk F,Verheij J,Li RJE,Schetters STT,Meijer LL,Le Large TYS,Driehuis E,Clevers H,Bruijns SCM,O'Toole T,van Vliet SJ,Bijlsma MF,Wuhrer M,Kazemier G,Giovannetti E, Garcia-Vallejo JJ,van Kooyk Y

**Year** 2022

Species Human

Dilution 1:1000

Journal of immunology (Baltimore, Md.: 1950)

Oligodeoxynucleotides differentially modulate activation of TLR7 and TLR8 by imidazoquinolines.

"AHC3419 was used in ELISA to investigate the cross-talk between oligodeoxynucleotides, immune response modifiers, TLR7, and TLR8"

Authors: Gorden KK,Qiu X,Battiste JJ,Wightman PP,Vasilakos JP,Alkan SS

**Year** 2006

Dilution 1:200

View more ELISA references on thermofisher.com

### More applications with references on thermofisher.com

For Research Use Only, Not for use in diagnostic procedures. Not for resale without express authorization, Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT.

BUYER'S EXCLUSIVE REMEDY FOR NON-CORNORMOR PRODUCTS DURING THE WARRANTY PERFIGIO IS LIMITAGE. ACCIDENT OR FOR PAIR. REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS, AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPLAR, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS, AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPLAR, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS, AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPLAR, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS OF THE PRODUCTS OF THE PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses