## CD27 Monoclonal Antibody (O323), NovaFluor™ Blue 660-120S, eBioscience™

## **Product Details**

Size	100 Tests
Host/Isotype	Mouse / IgG1, kappa
Class	Monoclonal
Туре	Antibody
Clone	O323
Conjugate	NovaFluor™ Blue 660-120S
Excitation/Emission Max	492/665 nm
Form	Liquid
Concentration	5 µL/Test
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_2896655

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	5 μL (0.8 μg)/test	-

## **Product Specific Information**

Description: The O323 monoclonal antibody reacts with human CD27, a lymphocyte-specific member of the TNFR superfamily. CD27 is expressed by a subset of thymocytes and virtually all mature T cells and is upregulated upon T-cell stimulation. CD27 binds to CD70, and through this interaction, plays an important role in T cell-B cell interaction.

Applications Reported: The O323 antibody has been reported for use in flow cytometric analysis.

Applications Tested: The O323 antibody has been pre-titrated and tested by flow cytometric analysis of normal human peripheral blood cells. This can be used at 5  $\mu$ L (0.8  $\mu$ g) per test. A test is defined as the amount ( $\mu$ g) of antibody that will stain a cell sample in a final volume of 100  $\mu$ L. Cell number should be determined empirically but can range from 10^5 to 10^8 cells /test.

Each NovaFluor conjugate or kit is shipped with CellBlox Blocking Buffer. Use this buffer whenever staining with NovaFluor conjugates, including single-color compensation controls using cells. Whenever possible, we recommend adding CellBlox Blocking Buffer to antibody cocktails/master mixes prior to combining with cells. Add 5  $\mu$ L per sample (regardless of the number of NovaFluors in your panel) to use the antibody cocktail as intended. For single-color controls, use 5  $\mu$ L of CellBlox Blocking Buffer per 100 $\mu$ L of cell sample containing 10^3 to 10^8 cells.

Excitation: 509 nm; Emission: 665 nm; Laser: 488 nm (Blue) Laser

NovaFluor conjugates are based on Phiton<sup>™</sup> technology utilizing novel nucleic acid dye structures that allow for engineered fluorescent signatures with consideration for spillover and spread impacts. Learn more

1

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production dcumentation, specifications and/or substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production dcumentation, the warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation dcumentation, the warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation of sumple material and product will conform to such model or sample (unrunised by properly trained individuals. Unless otherwise stated in the Documentation of the Documentation of the Documentation of the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample (unrunised to by properly trained individuals. No THER WARRANTES, EXPRESS OR INNERD, EXE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTES OF MERCHANTABILITY, FITNESS FOR ANY PARTICURS DENON INFRINCEMENT. BUYER'S EXCLUSIVE REMEMPT FOR NON-CONFORMING PRODUCTS SOLE OPTION. THERE IS NO BUILGATION TO REPARI, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS SOLE OPTION. THERE IS NO BUILGATION TO REPARI, REPLACE OR REFUND FOR THE NON-CONFORMING PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE RODUCTS. Unless otherwise expressly stated on the documentation or animals.