




Nestin Monoclonal Antibody (196908), APC

Product Details	
Size	100 Tests
Species Reactivity	Human
Published Species	Human
Host/Isotype	Mouse / IgG1
Class	Monoclonal
Туре	Antibody
Clone	196908
Conjugate	APC
Excitation/Emission Max	651/660 nm
Immunogen	NS0 mouse myeloma cell line transfected with human Nestin
Form	Liquid
Purification	Protein A/G
Storage buffer	Buffered saline solution with BSA
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_2608686

Applications	Tested Dilution	Publications
Flow Cytometry (Flow)	10 μL per million cells	3 Publications

Product Images For Nestin Monoclonal Antibody (196908), APC



Nestin Antibody (MA5-23650) in Flow

Flow cytometry of Nestin in NTera,2 human testicular embryonic carcinoma cell line. Samples were incubated in Nestin monoclonal antibody (Product # MA5-23650) or isotype control.

☐ 3 References

Flow Cytometry (3)

Cells

Combining HDAC and MEK Inhibitors with Radiation against Glioblastoma-Derived Spheres.

"Published figure using Nestin monoclonal antibody (Product # MA5-23650) in Flow Cytometry" Authors: Essien El,Hofer TP,Atkinson MJ,Anastasov N

Year 2022

Stem cell research

Generation of two H1 hESC sublines carrying deletions of RB1 exon 1 /promoter in heterozygous or compound heterozygous state.

"MA5-23650 was used in Flow cytometry/Cell sorting to modify hESCs, removing the RB1 exon 1 using the CRISPR /Cas9 nickase system."

Authors: Menges J, Cremanns M, Steenpass L

Year 2019

Species Human

View more Flow references on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"), No claim of suitability for use in applications regulated by EDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is imited to one year from date of subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample, NO OTHER WARRANTIES, EXPRESS OR IMPLED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, ETITIESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT.

BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMINO PRODUCTS DURING the WARRANTIES OF DEPAIR, REPLACE OR REFUND FOR THE NON-CONFORMINO PRODUCTS AS THE RESULT OF () ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (ii) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (iii) USE OF THE PRODUCTS IN A MANNER FOR WHIch THEY WERE NOT DESIGNED, OR (iv) IMPROPER'S TORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, or vivo or in vivo therapeutic uses, or any type of consumption to human or animal to human or