

# Phospho-p70 S6 Kinase (Thr421) Polyclonal Antibody


Product Details	
Size	100 µL
Species Reactivity	Human, Mouse, Rat
Published Species	Mouse
Host/Isotype	Rabbit / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	Peptide sequence around phosphorylation site of threonine 421 (P-R-T(p)-P-V) derived from Human p70 S6 Kinase.
Form	Liquid
Concentration	1 mg/mL
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.4, with 50% glycerol
Contains	0.02% sodium azide
Storage conditions	-20°C
RRID	AB_2554341

Applications	Tested Dilution	Publications
Western Blot (WB)	1:500-1:1,000	-
Immunohistochemistry (IHC)	-	1 Publication
Immunohistochemistry (Paraffin) (IHC (P))	1:50-1:100	1 Publication
Immunocytochemistry (ICC/IF)	1:100-1:200	-


## Product Specific Information

A suggested positive control for Western blot is Jurkat cells; suggested positive control for IHC is human breast carcinoma; suggested positive control for ICC/IF is MCF-7 cells.

Product Images For Phospho-p70 S6 Kinase (Thr421) Polyclonal Antibody



**Phospho-p70 S6 Kinase (Thr421) Antibody (PA5-37733) in ICC/IF**  
Immunocytochemical analysis of Phospho-p70 S6 Kinase (Thr421) in methanol-fixed MCF cells, using Phospho-p70 S6 Kinase (Thr421) Polyclonal Antibody (Product # PA5-37733).



**Phospho-p70 S6 Kinase (Thr421) Antibody (PA5-37733) in IHC (P)**  
Immunohistochemical analysis of Phospho-p70 S6 Kinase (Thr421) in paraffin-embedded Human breast carcinoma tissue using Phospho-p70 S6 Kinase (Thr421) Polyclonal Antibody (Product # PA5-37733) (left) or the same antibody preincubated with blocking peptide (right).

2 References

Immunohistochemistry (1)

eLife	Year 2020
Increased mTOR activity and metabolic efficiency in mouse and human cells containing the African-centric tumor-predisposing p53 variant Pro47Ser.	Species Mouse
"PA5-37733 was used in Immunohistochemistry (Paraffin) to show that mice containing the S47 variant display increased mTOR activity and oxidative metabolism, as well as larger size, improved metabolic efficiency, and signs of superior fitness."	Dilution 1:100
Authors: Gnanapradeepan K,Leu JI,Basu S,Barnoud T,Good M,Lee JV,Quinn WJ,Kung CP,Ahima R,Baur JA,Wellen KE,Liu Q,Schug ZT,George DL,Murphy ME	

Immunohistochemistry (Paraffin) (1)

eLife	Year 2020
Increased mTOR activity and metabolic efficiency in mouse and human cells containing the African-centric tumor-predisposing p53 variant Pro47Ser.	Species Mouse
"PA5-37733 was used in Immunohistochemistry (Paraffin) to show that mice containing the S47 variant display increased mTOR activity and oxidative metabolism, as well as larger size, improved metabolic efficiency, and signs of superior fitness."	Dilution 1:100
Authors: Gnanapradeepan K,Leu JI,Basu S,Barnoud T,Good M,Lee JV,Quinn WJ,Kung CP,Ahima R,Baur JA,Wellen KE,Liu Q,Schug ZT,George DL,Murphy ME	

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON-INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.