

VEGF Receptor 2 Polyclonal Antibody

Product Details	
Size	100 µg
Species Reactivity	Mouse
Published Species	Human, Mouse
Host/Isotope	Rabbit / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	Synthetic peptide derived from the N-terminal region of the mouse vascular endothelial growth factor receptor-2 (VEGFR-2, FLK-1, KDR) protein.
Form	Liquid
Concentration	0.25 mg/mL
Purification	Antigen affinity chromatography
Storage buffer	PBS, pH 7.4
Contains	0.1% sodium azide
Storage Conditions	-20°C
RRID	AB_2533237

Applications	Tested	Dilution	Published
Immunofluorescence (IF)	-		1 Publication
Miscellaneous PubMed (MISC)	-		1 Publication
Western Blot (WB)	✓	Assay Dependent	1 Publication
Immunohistochemistry (Frozen) (IHC (F))	-		1 Publication
ELISA (ELISA)	✓	Assay Dependent	

4 References

Immunofluorescence (1)

Journal of cellular physiology

Expression of receptors for VEGFs on normal human thyroid follicular cells and their role in follicle formation.

Authors: Susarla R,Gonzalez AM,Watkinson JC,Eggo MC

Species

Human
Not Applicable

Dilution

Not Cited
Not Cited

Year

2012

Miscellaneous PubMed (1)

Journal of cellular physiology

Expression of receptors for VEGFs on normal human thyroid follicular cells and their role in follicle formation.

Authors: Susarla R,Gonzalez AM,Watkinson JC,Eggo MC

Species

Human
Not Applicable

Dilution

Not Cited
Not Cited

Year

2012

Western Blot (1)

Journal of cellular physiology

Expression of receptors for VEGFs on normal human thyroid follicular cells and their role in follicle formation.

Authors: Susarla R,Gonzalez AM,Watkinson JC,Eggo MC

Species

Human
Not Applicable

Dilution

Not Cited
Not Cited

Year

2012

Immunohistochemistry (Frozen) (1)

Genes and development

Coronary vessel development requires activation of the TrkB neurotrophin receptor by the Wilms' tumor transcription factor Wt1.

Authors: Wagner N,Wagner KD,Theres H,Englert C,Schedl A,Scholz H

Species

Mouse

Dilution

Not Cited

Year

2005

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.