

CD4 Monoclonal Antibody (RM4-5), Alexa Fluor™ 700, eBioscience™

Product Details	
Size	100 µg
Species Reactivity	Mouse
Published Species	Mouse, Human
Host/Isotype	Rat / IgG2a, kappa
Recommended Isotype Control	Rat IgG2a kappa Isotype Control (eBR2a), Alexa Fluor™ 700, eBioscience™
Class	Monoclonal
Type	Antibody
Clone	RM4-5
Conjugate	Alexa Fluor™ 700
Excitation/Emission Max	696/719 nm
Form	Liquid
Concentration	0.2 mg/mL
Purification	Affinity chromatography
Storage buffer	PBS, pH 7.2
Contains	0.09% sodium azide
Storage conditions	4° C, store in dark, DO NOT FREEZE!
RRID	AB_494000

Applications	Tested Dilution	Publications
Immunohistochemistry (IHC)	-	10 Publications
Immunohistochemistry (PFA fixed) (IHC (PFA))	-	1 Publication
Immunohistochemistry (Frozen) (IHC (F))	-	1 Publication
Immunocytochemistry (ICC/IF)	-	4 Publications
Flow Cytometry (Flow)	0.125 µg/test	136 Publications

Product Specific Information

Description: The RM4-5 monoclonal antibody reacts with the mouse CD4 molecule, a 55 kDa cell surface receptor expressed by a majority of thymocytes, subpopulation of mature T cells and dendritic cells. CD4 binds to MHC class II on the surface of antigen presenting cells and plays an important role both in T cell development and in optimal functioning of mature T cells. In T cells, CD4 associates with protein tyrosine kinase p56lck through its cytoplasmic tail. Binding of RM4-5 is blocked by GK1.5.

Applications Reported: This RM4-5 antibody has been reported for use in flow cytometric analysis.

Applications Tested: This RM4-5 antibody has been tested by flow cytometric analysis of mouse thymocyte and splenocyte suspensions. This can be used at less than or equal to 0.125 µg per test. A test is defined as the amount (µg) of antibody that will stain a cell sample in a final volume of 100 µL. Cell number should be determined empirically but can range from 10⁵ to 10⁸ cells/test. It is recommended that the antibody be carefully titrated for optimal performance in the assay of interest.


Alexa Fluor® 700 emits at 723 nm and can be excited with the red laser (633 nm). Most instruments will require a 685 LP

mirror and 710/20 filter. Please make sure that your instrument is capable of detecting this fluorochrome.

Excitation: 633-647 nm; Emission: 723 nm; Laser: Red Laser.


Filtration: 0.2 µm post-manufacturing filtered.

Product Images For CD4 Monoclonal Antibody (RM4-5), Alexa Fluor™ 700, eBioscience™



CD4 Antibody (56-0042-82)

Staining of mouse splenocytes and bone marrow cells. Right: As expected based on known relative expression patterns, CD4 clone RM4-5 stains a subset of splenocytes and does not stain any bone marrow cells. Details: Balb/c bone marrow cells (left) and splenocytes (middle) were surface stained with CD4 (clone RM4-5) followed by staining with 7-AAD. Viable bone marrow cells in the lymphoid (blue histogram) and myeloid (purple histogram) gates and viable splenocytes (orange histogram) were used for analysis. {RE}



CD4 Antibody (56-0042-82) in Flow

Staining of mouse splenocytes with Anti-Mouse CD3 FITC (Product # 11-0032-82) and 0.06 µg of Rat IgG2a K Isotype Control Alexa Fluor® 700 (Product # 56-4321-80) (left) or 0.06 µg of Anti-Mouse CD4 Alexa Fluor® 700 (right). Cells in the lymphocyte gate were used for analysis.

[View more figures on thermofisher.com](#)

Immunohistochemistry (10)

Cells	Year 2023
Complete Freund's Adjuvant Induces a Fibroblast-like Synoviocytes (FLS) Metabolic and Migratory Phenotype in Resident Fibroblasts of the Inoculated Footpad at the Earliest Stage of Adjuvant-Induced Arthritis.	
"Published figure using CD4 monoclonal antibody (Product # 56-0042-82) in Immunohistochemistry"	
Authors: González-Chávez SA,Chaparro-Barrera E,Alvarado-Jáquez MF,Cuevas-Martínez R,Ochoa-Albíztegui RE,Pacheco-Tena C	
Frontiers in immunology	Year 2022
Black raspberry extract inhibits regulatory T-cell activity in a murine model of head and neck squamous cell carcinoma chemoprevention.	
"Published figure using CD4 monoclonal antibody (Product # 56-0042-82) in Immunohistochemistry"	
Authors: Ryan NM,Lamenza FF,Upadhaya P,Pracha H,Springer A,Swingler M,Siddiqui A,Oghumu S	

View more IHC references on thermofisher.com

Immunohistochemistry (PFA fixed) (1)

Immunity	Year 2019
T Cell Recruitment to the Intestinal Stem Cell Compartment Drives Immune-Mediated Intestinal Damage after Allogeneic Transplantation.	
"Published figure using CD4 monoclonal antibody (Product # 56-0042-82) in Immunocytochemistry"	
Authors: Fu YY,Egorova A,Sobieski C,Kuttiyara J,Calafiore M,Takashima S,Clevers H,Hanash AM	

Immunohistochemistry (Frozen) (1)

Biology	Year 2021
Increased Levels of IL-16 in the Central Nervous System during Neuroinflammation Are Associated with Infiltrating Immune Cells and Resident Glial Cells.	
"Published figure using CD4 monoclonal antibody (Product # 56-0042-82) in Immunohistochemistry (Frozen)"	
Authors: Hridi SU,Barbour M,Wilson C,Franssen AJ,Harte T,Bushell TJ,Jiang HR	

More applications with references on thermofisher.com

- ICC/IF (4)
- Flow (136)

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON-INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.