

CD34 Monoclonal Antibody (581), FITC

Product Details	
Size	500 µL
Species Reactivity	Human
Published Species	Human
Host/Isotype	Mouse / IgG1
Class	Monoclonal
Type	Antibody
Clone	581
Conjugate	FITC
Excitation/Emission Max	498/517 nm
Immunogen	Human CD34.
Form	Liquid
Concentration	0.2 mg/mL
Purification	purified
Storage buffer	PBS with 4mg/mL BSA
Contains	0.1% sodium azide
Storage conditions	4° C, store in dark
RRID	AB_2536499

Applications	Tested Dilution	Publications
Immunohistochemistry (IHC)	-	1 Publication
Flow Cytometry (Flow)	5 µL for 1x10^6 cells	13 Publications

Immunohistochemistry (1)

Pathology international Phyllodes tumor of the seminal vesicle: case report and literature review. "Published figure using CD34 monoclonal antibody (Product # CD34-581-01) in Immunohistochemistry" Authors: Son HJ,Jeong YJ,Kim JH,Chung MJ	Year 2004
--	----------------------

Flow Cytometry (13)

Experimental and therapeutic medicine Bcl-xL mutant promotes cartilage differentiation of BMSCs by upregulating TGF-/BMP expression levels. "Published figure using CD34 monoclonal antibody (Product # CD34-581-01) in Flow Cytometry" Authors: Xiao K,Yang L,Xie W,Gao X,Huang R,Xie M	Year 2021 Species Human
---	---

Stem cells international Impact of Tissue Harvesting Sites on the Cellular Behaviors of Adipose-Derived Stem Cells: Implication for Bone Tissue Engineering. "CD34-581-01 was used in Flow cytometry/Cell sorting to quantitatively compare the proliferation and osteogenic capabilities of adipose-derived stem cells from different harvesting sites." Authors: Rezai Rad M,Bohlooli M,Akhavan Rahn timer M,Anbarlou A,Nazeman P,Khojasteh A	Year 2020 Species Human
--	---

View more Flow references on thermofisher.com

More applications with references on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (i) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (ii) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (iii) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (iv) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.