

# C/EBP alpha Monoclonal Antibody (15C8)

Product Details	
Size	100 µg
Species Reactivity	Human, Mouse, Rat
Published Species	Non-human primate, Mouse, Human
Host/Isotype	Mouse / IgG1
Class	Monoclonal
Type	Antibody
Clone	15C8
Conjugate	Unconjugated
Immunogen	Synthetic peptide corresponding to residues M(1) E S A D F Y E V E P R P P(14) of mouse C/EBP alpha.
Form	Liquid
Concentration	1 mg/mL
Purification	Protein G
Storage buffer	PBS with 1mg/mL BSA
Contains	0.05% sodium azide
Storage conditions	-20° C, Avoid Freeze/Thaw Cycles
RRID	AB_2077903

Applications	Tested Dilution	Publications
Western Blot (WB)	2 µg/mL	5 Publications
Immunocytochemistry (ICC/IF)	-	1 Publication
Immunoprecipitation (IP)	Assay-dependent	2 Publications
ChIP assay (ChIP)	-	2 Publications

## Product Specific Information

MA1-825 detects C/EBP alpha protein in human, mouse and rat samples.

MA1-825 has successfully been used in Western blot and immunoprecipitation procedures. By Western blot, this antibody detects a 45 kDa protein representing C/EBP alpha in transfected hu293T cells. MA1-825 is specific to the amino terminal region of C/EBP alpha.

The MA1-825 immunogen is a synthetic peptide corresponding to residues M(1) E S A D F Y E V E P R P P(14) of mouse C /EBP alpha. This peptide (Cat. # PEP-242) is available for use in neutralization and control experiments.

Fig. 1

**C/EBP alpha Antibody (MA1-825) in WB**  
Western blot detection of C/EBP alpha on transfected 293T cells using Product # MA1-825.

C/EBP alpha



[View more figures on thermofisher.com](https://www.thermofisher.com)

Western Blot (5)

Scientific reports	Year 2021
Antibesity and antidiabetic effects of the dairy bacterium <i>Propionibacterium freudenreichii</i> MJ2 in high-fat diet-induced obese mice by modulating lipid metabolism.	Species Mouse
"Published figure using C/EBP alpha monoclonal antibody (Product # MA1-825) in Western Blot"	Dilution 1:500
Authors: An M,Park YH,Lim YH	
Adipocyte	Year 2020
The histone methyltransferase Suv39h regulates 3T3-L1 adipogenesis.	Species Mouse
"Published figure using C/EBP alpha monoclonal antibody (Product # MA1-825) in Western Blot"	
Authors: Jing J,Li F,Zha L,Yang X,Wu R,Wang S,Xue B,Shi H	

View more WB references on thermofisher.com

Immunocytochemistry (1)

The EMBO journal	Year 2010
Nucleolar retention of a translational C/EBPalpha isoform stimulates rDNA transcription and cell size.	Species Human Non-human primate
"MA1-825 was used in immunocytochemistry to investigate the role of C/EBP-alpha on the regulation of rRNA synthesis and cell growth and proliferation"	
Authors: Müller C,Bremer A,Schreiber S,Eichwald S,Calkhoven CF	

Immunoprecipitation (2)

The EMBO journal	Year 2010
Nucleolar retention of a translational C/EBPalpha isoform stimulates rDNA transcription and cell size.	Species Human Non-human primate
"MA1-825 was used in immunocytochemistry to investigate the role of C/EBP-alpha on the regulation of rRNA synthesis and cell growth and proliferation"	
Authors: Müller C,Bremer A,Schreiber S,Eichwald S,Calkhoven CF	

More applications with references on thermofisher.com

ChIP (2)

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON-INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.