

LOX Polyclonal Antibody

Product Details

Size	100 µL
Species Reactivity	Human, Mouse, Pig, Rat
Host/Isotype	Rabbit / IgG
Class	Polyclonal
Туре	Antibody
Conjugate	Unconjugated
Immunogen	Synthetic peptide between amino acids 200-300 of the human LOX protein.
Form	Liquid
Concentration	1 mg/mL
Purification	Antigen affinity chromatography
Storage buffer	PBS
Contains	0.02% sodium azide
Storage conditions	-20° C, Avoid Freeze/Thaw Cycles
RRID	AB_794775

Applications	Tested Dilution	Publications
Western Blot (WB)	0.5-4 μg/mL	4 Publications
Immunohistochemistry (Paraffin) (IHC (P))	1:100-1:250	-
Immunocytochemistry (ICC/IF)	1:100	-

Product Specific Information


PA1-16953 is specific for the ~58 kDa LOX protein.

This is a highly glycosylated protein and therefore, other bands may be seen at ~32 kDa (non-glycosylated form) 50 kDa and 58 kDa (glycosylated forms).

Suggested positive control: human and mouse kidney lysate, antigen standard for LOX (transient overexpression lysate), human kidney protein, mouse kidney protein.

1

Product Images For LOX Polyclonal Antibody




LOX Antibody (PA1-16953) in WB


Western blot was performed using Anti-LOX Rabbit Polyclonal Antibody (Product # PA1-16953) and 47 kDa (canonical) and 60 kDa (glycosylated) bands corresponding to isoforms of LOX were observed in tissues tested except Mouse Spleen which is reported to be negative. Tissue extracts (30 µg lysate) of Mouse Heart (Lane 1), Mouse Kidney (Lane 2) and Mouse Spleen (Lane 3) were electrophoresed using Novex® NuPAGE® 4-12% Bis-Tris Protein Gel (Product # NP0322BOX). Resolved proteins were then transferred onto a nitrocellulose membrane (Product # IB23001) by iBlot® 2 Dry Blotting System (Product # IB21001). The blot was probed with the primary antibody (2 µg/mL) and detected by chemiluminescence with Goat anti-Rabbit IgG (Heavy Chain) Superclonal™ Recombinant Secondary Antibody, HRP (Product # A27036, 1:4000 dilution) using the iBright FL 1000 (Product # A32752). Chemiluminescent detection was performed using Novex® ECL Chemiluminescent Substrate Reagent Kit (Product # WP20005).

LOX Antibody (PA1-16953)

Antibody specificity was demonstrated by detection of differential basal expression of the target across tissues tested owing to their inherent genetic constitution. Relative expression of LOX was observed in Mouse Heart (Lane 1) and Mouse Kidney (Lane 2) that express LOX protein, in comparison to, Mouse Spleen that doesnt express LOX protein using Anti-LOX Rabbit Polyclonal Antibody (Product # PA1-16953) in Western Blot. {RE}







LOX Antibody (PA1-16953) in IHC (P)

Immunohistochemical analysis of LOX in mouse stomach. Samples were incubated in LOX polyclonal antibody (Product # PA1-16953).

View more figures on thermofisher.com

□ 4 References

Western Blot (4)

Frontiers in pharmacology DuCLOX-2/5 Inhibition Attenuates Inflammatory Response and Induces Mitochondrial Apoptosis for Mammary Gland Chemoprevention.		
"Published figure using LOX polyclonal antibody (Product # PA1-16953) in Western Blot"		
Authors: Gautam S,Rawat AK,Sammi SR,Roy S,Singh M,Devi U,Yadav RK,Singh L,Rawat JK,Ansari MN,Saeedan AS, Kumar D,Pandey R,Kaithwas G		
Journal of cell science Defective signaling, osteoblastogenesis and bone remodeling in a mouse model of connexin 43 C-terminal truncation.	Year 2017	

"Published figure using LOX polyclonal antibody (Product # PA1-16953) in Western Blot" Authors: Moorer MC,Hebert C,Tomlinson RE,Iyer SR,Chason M,Stains JP

View more WB references on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Productin documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product till conform to such model or sample. No OTHER WARRANTES, EXPRESS OR IMPLED, ARRANTES DE MERCHANTABILITY, THYRESS FOR ANY PARTICULAR PURPRESS, OR NIMELON, ARRANTES DE MERCHANTABILITY, THYRESS FOR ANY PARTICULAR PURPRESS, OR NIMELON ARRANTY PERIOD EXPRESS OR REFUND FOR THE NON-CONFORMING PRODUCTS) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCTS () A SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCTS) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCTS () A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS AS THE RESULT OF () ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FLOR ON CONFORMING PRODUCTS) IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS SO THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE