

# Androgen Receptor Polyclonal Antibody

Product Details	
Size	500 µL
Species	Human, Mouse, Rat
Published Species	Rat, Human
Expression System	Rabbit / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	A synthetic peptide derived from N-terminus of rat AR
Form	Liquid
Concentration	1 mg/mL
Purification	Protein A
Storage buffer	PBS, pH 7.4, with 0.2% BSA
Contains	0.09% sodium azide
Storage Conditions	4° C
RRID	AB_10984711


Applications	Tested Dilution	Publications
Immunocytochemistry (ICC)	2 µg/mL	2 Publications
Immunofluorescence (IF)	2 µg/mL	-
Immunoprecipitation (IP)	10 µg/mL	-
Western Blot (WB)	5 µg/mL	-
Immunohistochemistry (IHC)	-	1 Publication

## Product Specific Information

PA5-16363 targets Androgen Receptor in IP and WB applications and shows reactivity with mouse, Rat, and Human samples.

The PA5-16363 immunogen is a synthetic peptide derived from N-terminus of rat AR.


## Advanced Verification Data



### Androgen Receptor Antibody (PA5-16363)


IP-MS enrichment of AR (LFQ intensity): AR was enriched 211-fold from LNCAP lysate compared to background proteins, using the optimized IP-MS workflow with Pierce MS-Compatible Magnetic IP Kit protein A/G (Product # 90409) and AR antibody (Product # PA5-16363). The STRING database ([www.string-db.org](http://www.string-db.org)) was used to identify the protein interactor list. See more information on IP-MS verification of antibody selectivity. IP-MS validation info.

## Product Images For Androgen Receptor Polyclonal Antibody



### Androgen Receptor Antibody (PA5-16363) in IF

Immunofluorescence analysis of Androgen receptor was performed using 70% confluent log phase PC-3 cells. The cells were fixed with 4% paraformaldehyde for 10 minutes, permeabilized with 0.1% Triton™ X-100 for 10 minutes, and blocked with 1% BSA for 1 hour at room temperature. The cells were labeled with Androgen Receptor Rabbit Polyclonal Antibody (Product # PA5-16363) at 2µg/mL in 0.1% BSA and incubated for 3 hours at room temperature and then labeled with Goat anti-Rabbit IgG (H+L) Superclonal™ Secondary Antibody, Alexa Fluor® 488 conjugate (Product # A27034) at a dilution of 1:2000 for 45 minutes at room temperature (Panel a: green). Nuclei (Panel b: blue) were stained with SlowFade® Gold Antifade Mountant with DAPI (Product # S36938). F-actin (Panel c: red) was stained with Alexa Fluor® 555 Rhodamine Phalloidin (Product # R415, 1:300). Panel d represents the merged image showing cytoplasmic localization. Panel e shows the no primary antibody control. The images were captured at 60X magnification.



### Androgen Receptor Antibody (PA5-16363) in IHC

Formalin-fixed, paraffin-embedded Human breast cancer tissue stained with Androgen Receptor using peroxidase-conjugate and AEC chromogen.

[View more figures on thermofisher.com](http://thermofisher.com)

## Immunocytochemistry (2)

### Scientific reports

#### Altered hormonal milieu and dysregulated protein expression can cause spermatogenic arrest in ectopic xenografted immature rat testis.

"PA5-16363 was used in Western Blotting to investigate the underlying causes of spermatogenic arrest in rats."

Authors: Goel S,Minami N

**Species**  
Rat

**Dilution**  
1:100

**Year**  
2019

### PloS one

#### Targeting androgen receptor/Src complex impairs the aggressive phenotype of human fibrosarcoma cells.

"PA5-16363 was used in immunocytochemistry to study the inhibition of the aggressive phenotype of a human fibrosarcoma cell line by disruption of the interaction between the AR and Src tyrosine kinase"

Authors: Castoria G,Giovannelli P,Di Donato M,Hayashi R,Arra C,Appella E,Auricchio F,Migliaccio A

**Species**  
Human

**Dilution**  
1:100

**Year**  
2014

## Immunohistochemistry (1)

### Andrologia

#### Testosterone administration to adult rats differentially modulates androgen and oestrogen receptor- expression in reproductive organs and pituitary.

"PA5-16363 was used in immunohistochemistry and western blot to study the effect of testosterone administration on androgen and oestrogen receptor-alpha expression in reproductive organs and pituitary of rats"

Authors: Kaushik MC,Misro MM,Sehgal N,Nandan D

**Species**  
Rat

**Dilution**  
1:100

**Year**  
2012

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.