

Phospho-c-Jun (Ser73) Polyclonal Antibody

Product Details	
Size	100 µL
Species Reactivity	Human, Mouse, Non-human primate, Rat
Host/Isotype	Rabbit / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	Synthetic phosphopeptide corresponding to residues around pSer73 of human c-Jun
Form	Liquid
Concentration	8 µg/mL
Purification	Antigen affinity chromatography
Storage buffer	0.01M HEPES, pH 7.5, with 0.15M NaCl, 100µg/mL BSA, 50% glycerol
Contains	no preservative
Storage conditions	-20°C
RRID	AB_10979275

Applications	Tested Dilution	Publications
Western Blot (WB)	1:1,000	-
Immunocytochemistry (ICC/IF)	1:100	-


Product Specific Information

It is not recommended to aliquot this antibody.

Product Images For Phospho-c-Jun (Ser73) Polyclonal Antibody


Phospho-c-Jun (Ser73) Antibody (PA5-17879) in WB

Western blot analysis was performed on whole cell extracts (30 µg lysate) of NIH/3T3 (Lane 1), NIH/3T3 treated with Anisomycin (10 µg/mL for 30 min) (Lane 2), NIH/3T3 treated with EGF (100 ng/mL for 15 min) (Lane 3), NIH/3T3 treated with TNF alpha (10 ng/mL for 15 min) (Lane 4), A-431 (Lane 5), A-431 treated with EGF (20 ng/mL for 60 min) (Lane 6) and A-431 treated with Afatinib (0.5 µM for 6 h) followed by EGF (20 ng/mL for 60 min) (Lane 7). The blot was probed with Anti-Phospho-c-Jun (Ser73) Polyclonal Antibody (Product # PA5-17819, 1:500 dilution) and detected by chemiluminescence using Goat anti-Rabbit IgG (Heavy Chain) Superclonal™ Secondary Antibody, HRP conjugate (Product # A27036, 0.25 µg/mL, 1:4,000 dilution). A 42 kDa band corresponding to Phospho-c-Jun (Ser63) was observed across the cell lines tested upon Anisomycin, EGF and TNF alpha treatment and was also observed to downregulate upon Afatinib treatment.




Phospho-c-Jun (Ser73) Antibody (PA5-17879)

Altered expression of target protein upon cell treatment demonstrates antibody specificity. Western blot analysis of Phospho-c-Jun (Ser73) Polyclonal Antibody (Product # PA5-17819) shows increased expression of proteins phosphorylated at the Serine residue in NIH/3T3 and A-431 cell lines upon Anisomycin, EGF and TNF alpha treatments. Phosphorylation is observed to be downregulated in A-431 cells upon treatment with Afatinib. {TM}



Phospho-c-Jun (Ser73) Antibody (PA5-17879) in ICC/IF

Immunofluorescent analysis of Phospho-c-Jun pSer73 in C6 cells, anisomycin treated, using a Phospho-c-Jun pSer73 polyclonal antibody (Product # PA5-17879) (green) and a beta-Tubulin monoclonal antibody (red). DNA is labeled using a fluorescent blue dye.



View more figures on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.