

LBP Polyclonal Antibody

Product Details	
Size	100 µL
Species Reactivity	Human
Host/Isotype	Rabbit / IgG
Class	Polyclonal
Type	Antibody
Conjugate	Unconjugated
Immunogen	Recombinant protein fragment containing a sequence corresponding to a region within amino acids 166 and 467 (P18428) of LBP
Form	Liquid
Concentration	0.46 mg/mL
Purification	Antigen affinity chromatography
Storage buffer	0.1M tris glycine, pH 7, with 10% glycerol
Contains	0.01% thimerosal
Storage conditions	Store at 4°C short term. For long term storage, store at -20°C, avoiding freeze/thaw cycles.
RRID	AB_11152918


Applications	Tested Dilution	Publications
Western Blot (WB)	1:500-1:3,000	1 Publication
Immunohistochemistry (Paraffin) (IHC (P))	1:100-1:1,000	-
Immunocytochemistry (ICC/IF)	-	1 Publication

Product Specific Information

Recommended positive controls: A431, H1299.


Store product as a concentrated solution. Centrifuge briefly prior to opening the vial.

Product Images For LBP Polyclonal Antibody



LBP Antibody (PA5-21642) in IHC (P)

Immunohistochemical analysis of paraffin-embedded human gastric N+T, using LBP (Product # PA5-21642) antibody at 1:100 dilution. Antigen Retrieval: EDTA based buffer, pH 8.0, 15 min.



LBP Antibody (PA5-21642) in WB

Western Blot using LBP Polyclonal Antibody (Product # PA5-21642). Sample (30 µg whole cell lysate). A: H1299. 7.5% SDS PAGE. LBP Polyclonal Antibody (Product # PA5-21642) diluted at 1:500.

[View more figures on thermofisher.com](https://www.thermofisher.com)

Western Blot (1)

International immunology	Year 2022
Identification of CD14 and lipopolysaccharide-binding protein as novel biomarkers for sarcoidosis using proteomics of serum extracellular vesicles.	
"Published figure using LBP polyclonal antibody (Product # PA5-21642) in Western Blot"	
Authors: Futami Y, Takeda Y, Koba T, Narumi R, Nojima Y, Ito M, Nakayama M, Ishida M, Yoshimura H, Naito Y, Fukushima K, Takimoto T, Edahiro R, Matsuki T, Nojima S, Hirata H, Koyama S, Iwahori K, Nagatomo I, Shirai Y, Suga Y, Satoh S, Futami S, Miyake K, Shiroyama T, Inoue Y, Adachi J, Tomonaga T, Ueda K, Kumanogoh A	

Immunocytochemistry (1)

Materials (Basel, Switzerland)	Year 2019
Adhesion Behaviour of Primary Human Osteoblasts and Fibroblasts on Polyether Ether Ketone Compared with Titanium under In Vitro Lipopolysaccharide Incubation.	
"Published figure using LBP polyclonal antibody (Product # PA5-21642) in Immunocytochemistry"	
Authors: Benz K, Schöbel A, Dietz M, Maurer P, Jackowski J	

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.