

Goat anti-Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody, BODIPY FL


Product Details	
Size	1 mg
Species Reactivity	Mouse
Host/Isotype	Goat / IgG
Class	Polyclonal
Type	Secondary Antibody
Conjugate	BODIPY® FL
Immunogen	Gamma Immunoglobins Heavy and Light chains
Form	Lyophilized
Purification	purified
Storage buffer	PBS, pH 7.5, with 1.5% BSA
Contains	0.01% thimerosal
Storage conditions	-20° C, store in dark
RRID	AB_2536429

Applications	Tested Dilution	Publications
Immunocytochemistry (ICC/IF)	1 µg/mL	-
Miscellaneous PubMed (Misc)	-	21 Publications


Product Specific Information

Storage and reconstitution: reconstitute using 0.5 mL deionized water. The product can be stored up to 2 weeks at 2-8°C. For longer storage, divide into single-use aliquots and freeze at or below -20°C. Frozen aliquots are stable for at least 6 months.


Product Images For Goat anti-Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody, BODIPY FL



Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody (B-2752) in ICC/IF
Immunofluorescence analysis of Goat anti-Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody, BODIPY FL (B2752) was performed using HeLa cells stained with alpha Tubulin (236-10501) Mouse Monoclonal Antibody (A11126). The cells were fixed with 4% paraformaldehyde for 10 minutes, permeabilized with 0.1% Triton™ X-100 for 10 minutes, blocked with 1% BSA for 1 hour and labeled with 2 µg/mL of mouse primary antibody for 3 hours at room temperature. Goat anti-Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody, BODIPY FL (B2752) was used at concentration of 1 µg/mL in phosphate buffered saline containing 0.2 % BSA for 45 minutes at room temperature, for detection of alpha Tubulin in the cytoplasm (Panel a: green). Nuclei (Panel b: blue) were stained with DAPI in SlowFade® Gold Antifade Mountant (S36938). F-actin was stained with Rhodamine Phalloidin (Product # R415, 1:300) (Panel c: red). Panel d represents the composite image. No nonspecific staining was observed with the secondary antibody alone (panel f), or with an isotype control (panel e). The images were captured at 60X magnification..



Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody (B-2752) in ICC/IF
FluoCells prepared slide #2 (Product # F-14781), which shows bovine pulmonary artery endothelial cells (BPAEC) that have been stained with an anti-tubulin mouse monoclonal antibody in conjunction with BODIPY FL goat anti-mouse IgG (Product # B-2752) for labeling microtubules, Texas Red-X phalloidin (Product # T7471) for labeling F-actin and DAPI (Product # D1306, Product # D3571, Product # D21490) for labeling nuclei. This multiple-exposure image was acquired using bandpass optical filter sets appropriate for DAPI, FITC and Texas Red dye.



Mouse IgG (H+L) Cross-Adsorbed Secondary Antibody (B-2752) in ICC/IF
FluoCells® prepared slide #2 (Product # F-14781) showing bovine pulmonary artery endothelial cells probed with anti-bovine alpha-tubulin mouse monoclonal 236-10501 (Product # A11126) and visualized with BODIPY® FL goat anti-mouse IgG (Product # B-2752). The actin filaments were then labeled with Texas Red®-X phalloidin (Product # T7471) and the nuclei were counterstained with DAPI (Product # D1306, D3571). The multiple-exposure image was acquired using bandpass filters appropriate for Texas Red® dye, FITC and DAPI.

[View more figures on thermofisher.com](https://www.thermofisher.com)

21 References

Miscellaneous PubMed (21)

The Journal of biological chemistry

A pool of extramitochondrial frataxin that promotes cell survival.

Authors: Condò I, Ventura N, Malisan F, Tomassini B, Testi R

Species
Not Applicable

Dilution
Not Cited

Year
2006

The Journal of cell biology

Shear stress-induced endothelial cell polarization is mediated by Rho and Rac but not Cdc42 or PI 3-kinases.

Authors: Wojciak-Stothard B, Ridley AJ

Species
Not Applicable

Dilution
Not Cited

Year
2003

[View more Misc references on thermofisher.com](#)

More applications with references on thermofisher.com

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization. Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Production documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample. NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS. Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to human or animals.