



LPAR5 Polyclonal Antibody

Catalog Number PA5-33899

Product data sheet

Details		Species Reactivity	
Size	50 µg	Species reactivity	Human
Host/Isotope	Rabbit / IgG	Tested Applications	
Class	Polyclonal	Immunohistochemistry (Paraffin) (IHC (P))	Dilution * 1-5 µg/mL
Type	Antibody	* Suggested working dilutions are given as a guide only. It is recommended that the user titrate the product for use in their own experiment using appropriate negative and positive controls.	
Immunogen	Synthetic 19 amino acid peptide from 3rd cytoplasmic domain of human LPAR5 / GPR92.		
Conjugate	Unconjugated		
Form	Liquid		
Concentration	1 mg/mL		
Storage Conditions	Maintain refrigerated at 2-8°C for up to 1 month. For long term storage store at -20°C		

Product specific information

Percent identity with other species by BLAST analysis: Human, Gorilla, Gibbon, Mouse, Panda (100%) Marmoset, Elephant, Bat, Bovine (95%).

Background/Target Information

Receptor for lysophosphatidic acid (LPA), a mediator of diverse cellular activities.

For Research Use Only. Not for use in diagnostic procedures. Not for resale without express authorization.